



Tissue Freezing

Freezing and Storage of Tissues

MACO BIOTECH Freezing – Storage Pots - REFERENCES CDR####AU

CE 0120



Medfor product

CE 0459



Former Maco Pharma product

Dear Customers,

Following the marketing cessation of our range of pots in November 2017 (references CDR1000AU, CDR1005AU, CDR2005AU, CDR7000AU, CDR8000AU) dedicated to the freezing and storage of tissues down to -80°C , Maco Pharma is happy to announce the establishment of a **new partnership with Medfor (UK company)**.

The references we will be able to distribute are described in the table overleaf.

This partnership will be effective from **February 2019** and will concern the following countries: Austria, Croatia, Czech Republic, Denmark, France, Germany, Iceland, Italy, Poland, Sweden.

For any additional information, please contact your local sales representative.

Best Regards,

The Biotherapy Marketing team

XA5IB01B 02/2019

Product Information





Tissue Freezing

	Reference	Raw materials	Bagging	Labelling	Volume (mL)		Sterilisation	Packaging	
					Int.	Ext.			
Maco Pharma Equivalents	CDR1000AU	SKA0050Z-K	Transparent polypropylene jars with white HDPE*tamper-evident screw caps	Double	2 pots labelled	50	100	Ethylene oxide	10 packs per case
		S241-K	Transparent polypropylene jars with white HDPE tamper-evident screw caps	Double	Outer pot labelled	2 x 50	2 x 100	Sterile R	30 packs per case
		S238-K	Transparent polypropylene jars with white HDPE tamper-evident screw caps	Double	Outer pot labelled	50	100	Sterile R	30 packs per case
	CDR1005AU CDR2005AU	SPA0050Z-K5	5 transparent polypropylene jars with white HDPE tamper-evident screw caps	Double	5 pots labelled	5 x 50		Ethylene oxide	10 packs per case
		SPA0100Z-K5	5 transparent polypropylene jars with white HDPE tamper-evident screw caps	Double	5 pots labelled	5 x 100		Ethylene oxide	10 packs per case
	CDR7000AU	SKA0300Z-K	Transparent polypropylene jars with white HDPE tamper-evident screw caps	Double	2 pots labelled	300	500	Ethylene oxide	10 packs per case
		S400-K	White polypropylene jars with white HDPE tamper-evident screw caps	Double	Outer pot labelled (pre-printed)				
		S401-K	White polypropylene jars with white HDPE tamper-evident screw caps	Double	No label				
		S411-K	White polypropylene jars with white HDPE tamper-evident screw caps	Double	Outer pot labelled (blank)				
	CDR8000AU	SKA0500Z-K	Transparent polypropylene jars with white HDPE tamper-evident screw caps	Double	2 pots labelled	500	1000	Ethylene oxide	8 packs per case
		S296-K	White polypropylene jars with white polypropylene screw caps	Double	Outer pot labelled (pre-printed)	500	750	Sterile R	10 packs per case
		S610-K	Transparent polypropylene jars with white HDPE tamper-evident screw caps	Double	Outer pot labelled (blank)	500	1000	Ethylene oxide	8 packs per case
	No Maco Pharma equivalent	S403-K	White polypropylene jars with white HDPE tamper-evident screw caps	Double	No label	100	300	Ethylene oxide	10 packs per case
		S402-K	White polypropylene jar with white HDPE tamper-evident screw cap, outer jar with white polypropylene screw cap	Double	Outer pot labelled (blank)	1250	2000	Sterile R	6 packs per case
		SKA0100Z-K	Transparent polypropylene jars with white HDPE tamper-evident screw caps	Double	2 pots labelled	100	300	Ethylene oxide	10 packs per case

*HDPE= High Density Polyethylene

The references **S241-K** and **S238-K** are not aimed to be frozen.

