



## evercare® Essentia Standard OP-Gowns, sterile

REF: 2522-01, 2523-01, 2524-01, 2525-01, 2527-01, 2528-01, 2529-01, 2530-01

evercare® Essentia OP-Gown is worn by a member of the surgical team to prevent transfer of infective agents. It is suitable for shorter procedures with small quantities of fluids, made in a soft SMS antistatic treated, with natural fluid repellent protection. All evercare® OP-Gowns complies with the requirements of EN 13795-1, provide an effective barrier to body fluids, skin particles and microbes and are designed for optimal comfort and ease of movement. The adjustable Velcro in the neck enables optimal fit while securing adjustments. High quality cuffs for improved comfort and tight fit round wrists. The gown is sterilized with EO-gas and packed in easy-to-open unit pouches. State-of-the-art packaging system with handy functional dispenser and labeling as per EN15223-1 and EN1041. Unit packs are supplied with barcoded removable labels for improved traceability.

### Key Features

- ▶ High breathability & high abrasion resistance
- ▶ SMS antistatic treated, with natural fluid repellent protection
- ▶ Ultrasonic heat sealed seams for sleeves and shoulders to increase comfort by reducing excessive materials causing discomfort
- ▶ Gown with overlapping back and adjustable Velcro in the neck enables optimal fit for individual adjustments
- ▶ Clear stamp with size and performance type
- ▶ Belt card with clear indications of sterile / non-sterile for secure donning Book-folding enabling unpacking and donning without risk for contamination SMS material with low linting



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2522-01	Essentia Standard OP-Gown, sterile, Blue	M	20 / 40
2523-01	Essentia Standard OP-Gown, sterile, Blue	L	20 / 40
2524-01	Essentia Standard OP-Gown, sterile, Blue	XL	18 / 36
2525-01	Essentia Standard OP-Gown, sterile, Blue	XXL	18 / 36
2527-01	Essentia Standard OP-Gown with hand towels, sterile, Blue	M	20 / 40
2528-01	Essentia Standard OP-Gown with hand towels, sterile, Blue	L	20 / 40
2529-01	Essentia Standard OP-Gown with hand towels, sterile, Blue	XL	18 / 36
2530-01	Essentia Standard OP-Gown with hand towels, sterile, Blue	XXL	18 / 36

Sales Unit: Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-05 EN ISO 10993-07 EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class Is
<b>Compliance with product standards:</b>	EN ISO 11135, EN 13795-1 Standard performance
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN 1041, EN ISO 11607-1

## Material

Product Part	Raw Material & Characteristics
<b>Body, sleeve, belt</b>	PP SMS nonwoven, blue
<b>Cuff</b>	PES, rib knitted fabric
<b>Neck ribbon</b>	PP Spunbond nonwoven, white
<b>Velcro band</b>	PA woven ribbon with hook or loop
<b>Belt card</b>	Cellulose pulp (ECF), white
<b>Wrapping crepe</b>	Pulp (ECF), white
<b>Sewing thread</b>	PES, white
<b>Hand towel</b>	Scrim Reinforced Tissue, white

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 Years

## Packing information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	Medical grade paper, PP/PE film	242 x 334
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	381 x 279 x 320
<b>Transport Carton</b>	Carton, TCF or ECF quality	574 x 401 x 330
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
2522-01	6438129120417	6438129320411	6438129520415
2523-01	6438129120424	6438129320428	6438129520422
2524-01	6438129120431	6438129320435	6438129520439
2525-01	6438129120448	6438129320442	6438129520446
2527-01	643812912046	643812932046	643812952046
2528-01	643812912047	643812932047	643812952047
2529-01	643812912048	643812932048	643812952048
2530-01	6438129120493	6438129320497	6438129520491

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.