


Intrauterine aspiration



INTRACURETTE SR

Curetta semirigida tipo "Karman"
per aspirazione intrauterina

Flexible intrauterine suction curette - "Karman" style

Cod.	Descrizione - Description	Colore Colour
720156	lung. 24 cm - length 24 cm Ø esterno 4 mm - O.D. 4 mm	bianco white
720157	lung. 24 cm - length 24 cm Ø esterno 5 mm - O.D. 5 mm	beige beige
720114	lung. 24 cm - length 24 cm Ø esterno 6 mm - O.D. 6 mm	grigio grey
720115	lung. 24 cm - length 24 cm Ø esterno 7 mm - O.D. 7 mm	giallo yellow
720116	lung. 24 cm - length 24 cm Ø esterno 8 mm - O.D. 8 mm	arancio orange
720117	lung. 24 cm - length 24 cm Ø esterno 9 mm - O.D. 9 mm	rosso red
720118	lung. 24 cm - length 24 cm Ø esterno 10 mm - O.D. 10 mm	verde green
720119	lung. 24 cm - length 24 cm Ø esterno 11 mm - O.D. 11 mm	azzurro cyan
720120	lung. 24 cm - length 24 cm Ø esterno 12 mm - O.D. 12 mm	blu blue

🇮🇹 Per interruzioni volontarie di gravidanza o a seguito di aborto spontaneo. La flessibilità del materiale permette un curettaggio delicato dopo l'evacuazione. Disponibile in 9 diverse misure con codice colore per una facile identificazione. Si connette direttamente al tubo di aspirazione.

🇺🇸 For uterine evacuation and curettage in voluntary pregnancy termination or following miscarriage. Flexibility of the material allows a gentle curettage after evacuation. It is available in 9 different colour-coded sizes for an easy identification. It can be connected directly to the suction tubing.

Scatola da 50 pezzi.
Fornito in confezione singola, monouso, sterilizzato a EtO e marcato CE.

50 pcs. per box.
Supplied in single pouch, single use, EtO sterilized and CE marked.

