

ETHICON

a Johnson & Johnson company



ETHILON*
(Polyamide 6) Suture

ETHILON*

Monofilament Polyamide 6

SUTURE

Description: **Monofilament** For smooth passage through tissue.
Synthetic For minimal tissue reaction.
Non-absorbable Provides prolonged tensile strength retention in tissue.

Color: Blue, black or clear.

Range: 11-0 to 2 (USP). Supplied as needed sutures armed with, for example, MULTIPASS* and PRIME* needles. ETHILON* is also available as ligatures.

Indications: ETHILON* (black and clear) sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, microsurgical and neurosurgical procedures. Blue ETHILON* sutures are intended for use in skin closure. Typical areas of use include:












- Cuticular sutures
- Nerve adaptation
- Ophthalmology

Full details in the Instruction for Use included in every package.



ETHILON* Monofilament Polyamide 6

3/8 Circle Reverse Cutting Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 C-2	13 mm	6-0	BLACK	45	250 microns	36	667H
		5-0	BLACK	45	250 microns	36	668H
 M-2	13 mm	5-0	BLACK	45	360 microns	36	644H
 FS-3	16 mm	6-0	BLACK	45		36	660H
		6-0	BLACK	45	VISI-BLACK*	36	EH7177H
		5-0	BLACK	45		36	1628H
		5-0	BLUE	45		36	EH7823H
		4-0	BLUE	45		36	EH7824H
 FS-2	19 mm	5-0	BLUE	45		36	EH7790H
		4-0	BLACK	45		36	662H
		4-0	BLUE	75		36	EH7143H
		3-0	BLUE	45		36	EH7792H
		2-0	BLUE	90		36	F3217
 FS-2S	19 mm	4-0	BLUE	45	Slimmed	36	EH7791H
 FS-1	24 mm	5-0	BLACK	45		36	EH7640H
		4-0	BLACK	45		36	1629H
		4-0	BLUE	75		36	FJ5272H
		3-0	BLACK	45		36	663H
		3-0	BLUE	75		36	EH7795H
		3-0	BLUE	90		36	F3221
		2-0	BLUE	45		36	EH7826H
		2-0	BLUE	75		36	EH7827H
 FS	26 mm	4-0	BLUE	45		36	FK5242H
		3-0	BLUE	45		36	ES163H
		2-0	BLUE	45		36	EH7796H
 FSL	30 mm	4-0	BLUE	75		36	1670H
		3-0	BLUE	45		36	EH7936H
		3-0	BLUE	75		36	EH7939H
		3-0	BLACK	45		36	S2455
		2-0	BLUE	75		36	EH7798H
 FSLX	36 mm	3-0	BLUE	75		36	1673H
		2-0	BLUE	75		36	1674H
		0	BLUE	75		36	EH7809H
 FSLY	45 mm	2-0	BLACK	100		12	W736G
 LR	76 mm	2	BLACK	3x50		12	3460TG






ETHILON* Monofilament Polyamide 6

3/8 Circle Reverse Cutting PRIME* Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
P-6	8 mm	7-0	BLACK	45		36	1647H
P-1	11 mm	7-0	BLACK	45		36	1696H
		7-0	BLACK	45	MULTIPASS*	36	MPE1696H
		6-0	BLACK	45		36	697H
		6-0	BLACK	45	MULTIPASS*	36	MPE697H
		5-0	BLACK	45		24	W1611T
P-3	13 mm	6-0	BLACK	45	MULTIPASS*	36	1698H
		5-0	BLACK	45	MULTIPASS*	36	698H
		4-0	UNDYED	45		36	691H
		4-0	BLACK	45		36	699H
		4-0	BLUE	45		36	EH7801H
		4-0	BLUE	45	MULTIPASS*	36	MPE7801H
		4-0	BLACK	45	MULTIPASS*	12	1665G
PS-3	16 mm	5-0	BLUE	45		36	EH7930H
		5-0	BLUE	45	MULTIPASS*	36	MPE7930H
		4-0	BLACK	45	MULTIPASS*	36	EH7761H
		4-0	BLUE	75		36	JD5273H
		5-0	BLUE	45		36	EH7804H
PS-2	19 mm	5-0	BLUE	45	MULTIPASS*	36	MPE7804H
		4-0	BLUE	45		36	EH7162H
		4-0	BLUE	45		24	W1620T
		4-0	BLACK	45	MULTIPASS*	36	MPE1667H
		4-0	BLUE	45	MULTIPASS*	36	MPE1855H
		4-0	UNDYED	45		12	1611G
		3-0	BLUE	45		24	W1621T
		3-0	BLUE	45	MULTIPASS*	36	MPE1621H
		4-0	BLACK	45		36	1667H
PS-1	24 mm	4-0	BLUE	45		36	EH7957H
		4-0	BLUE	75		36	EH7956H
		4-0	BLUE	75	MULTIPASS*	36	MPE7956H
		3-0	BLUE	45		24	W1857T
		3-0	BLUE	75		36	EH7954H
		3-0	BLUE	75	MULTIPASS*	36	MPE7954H
		3-0	BLUE	45	MULTIPASS*	36	MPE1625H
PS	26 mm	3-0	BLUE	45		24	W1625T
		3-0	BLUE	75		24	W1685T
		2-0	BLUE	75		24	W1626T
		3-0	BLUE	75		36	EH7113H
PSL	30 mm	3-0	BLUE	75		36	EH7113H

ETHILON* Monofilament Polyamide 6


3/8 Circle Conventional Cutting PRIME* Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PC-1	13 mm	6-0	BLACK	45		36	JM5140H
 PC-3	16 mm	6-0	BLACK	45	MULTIPASS*	12	1866G
		5-0	BLUE	45	MULTIPASS*	36	MPE1616H
		5-0	BLUE	45		24	W1616T
		5-0	BLACK	45		12	1865G
 PC-5	19 mm	4-0	BLUE	45	MULTIPASS*	36	MPE1856H
		4-0	BLUE	45		24	W1856T
 KP-3	90 mm	1	BLACK	2x50	With retention tubing	12	W2797
 KP-3 KP-3	90 mm	1	BLACK	75		12	W689

1/2 Circle Reverse Cutting Needles



Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 M-1	13 mm	5-0	BLACK	45		36	648H
 J-1	17 mm	5-0	BLACK	75		36	F2423H
 X-1	22 mm	3-0	BLUE	45		36	EH7943H
 CPX	48 mm	1	BLACK	100		12	W738

Straight Cutting Needles














Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 KS	60 mm	3-0	BLACK	75		36	627H
		2-0	BLACK	75		36	628H

ETHILON* Monofilament Polyamide 6

1/2 Circle Taper Point Needles



Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 MH-1	31 mm	0	BLACK	150	Loop	12	W739
 CTX	48 mm	1	BLACK	150	Loop	12	W748

3/8 Circle Taper Point Needles




Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 BV130-3	3,7 mm, 150 microns	9-0	BLACK	13		12	FG2819
 BV75-3	3,8 mm, 75 microns	10-0	BLACK	13		12	W2870
		9-0	BLACK	13		12	W2871
 BV50-3	3,8 mm, 50 microns	11-0	BLACK	13		12	W2881
 BV130-4	4,7 mm, 150 microns	9-0	BLACK	13		12	W2813
		8-0	BLACK	13		12	W2812
 BVH100-3	5 mm, 100 microns	10-0	BLACK	13	1/2 circle	12	W2860
 BV75-4	5,1 mm, 75 microns	10-0	BLACK	13		12	W2850
 BV-4	5,1 mm, 100 microns	10-0	BLACK	15		12	2814G
 BV100-4	5,1 mm, 100 microns	10-0	BLACK	13		12	W2830
		9-0	BLACK	13		12	W2829
		9-0	BLACK	15		12	EH7448G
 BV-2	6,5 mm, 150 microns	10-0	BLACK	15		12	2810G
		9-0	BLACK	15		12	2809G
		8-0	BLACK	15		12	2808G
 BV130-5	6,5 mm, 130 microns	10-0	BLACK	13		12	W2810
		8-0	BLACK	13		12	W2808
 BV130-5	6,5 mm, 130 microns	10-0	BLACK	30		12	W2801
 BV-1	9,3 mm, 250 microns	8-0	BLACK	45		12	W8170
 BV-1	9,3 mm, 250 microns	7-0	BLACK	45		36	EH7446H

ETHILON* Monofilament Polyamide 6


3/8 Circle CS ULTIMA* Spatula

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-140-6 	6,5 mm, 150 microns	10-0	BLACK	15		12	7060G
		10-0	BLACK	30		12	U7061
CS-140-6 	6,5 mm, 150 microns	10-0	BLACK	30		12	U7003



3/8 Circle MICROPOINT* Spatula

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-35-C 	6 mm, 150 microns	10-0	BLACK	15		12	EH7980G
GS-9 	6,5 mm, 200 microns	10-0	BLACK	30		12	1718G
		9-0	BLACK	30		12	1717G
TG140-6  TG140-6	6,5 mm, 150 microns	11-0	BLACK	30		12	W1780
		10-0	BLACK	30		12	W1770
		9-0	BLACK	30		12	W1769


1/2 Circle MICROPOINT* Reverse Cutting Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
G-7 	8 mm	8-0	BLACK	45		12	W1765


1/2 Circle MICROPOINT* Spatula

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
TG160-6 	5,5 mm, 150 microns	10-0	BLACK	30		12	W1756
TG175-8 	7 mm, 200 microns	8-0	BLACK	30		12	1716G

1/2 Circle CS ULTIMA* Spatula


Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-160-6 	5,5 mm, 150 microns	10-0	BLACK	30		12	U7000

1/4 Circle MICROPOINT* Spatula



Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
TG100-8 	6,5 mm, 220 microns	8-0	BLACK	30		12	B7714

ETHILON* Monofilament Polyamide 6

1/4 Circle Spatula Needles

	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
S-14 	8 mm	5-0	BLACK	20	36	MN5132H

1/4 Circle TAPERCUT* Needles

	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
V100-3 	3,5 mm, 100 microns	10-0	BLACK	13	12	W2899
		9-0	BLACK	13	12	W2898
V150-4 	4,5 mm, 150 microns	10-0	BLACK	15	12	W1701

Ligatures

	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
LIGAPAK*	2-0	BLACK	250	Reel	20	F4354