



Product description

Order number	05.1729.001
Product description	S-Monovette® 2.7 ml, w/o preparation, 66x11 mm, white EU/US code, paper label, 50/inner box sterile

Product characteristics

Type of collection	Venous
Type of preparation	Without preparation
Colour code	EU
Label/ Print	With paper label
Colour of print/label	Black
Cap	Membrane screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	Room temperature
Centrifugation time	10 min
Connection	Membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	2,7 ml
Diameter	11 mm
Length including cap	82 mm
Length excluding cap	66 mm
Length including cap and plunger	103 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	Transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	White
Membrane material	Rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product in accordance with MDD or IVDD	IVD
CE certified	CE
Sterilisation	sterile
Lot or serial number	batched

Packaging

Minimum order qty.	500
Type of packaging	Box
Pcs. / inner box	50
Pcs. / box	500
Pcs. / pallet	56000
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0,0126 cbm
Weight of case	2,45 KG

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





EAN of inner box	4038917133371
EAN of outer box	4038917024853

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

