

Technical Specifications

Brand:

**ZILLION
BLACK**

Description: Tissue forceps 1x2 teeth

Final dimensions:	14.5 cm
Composition:	High performance polymer.
Latex:	Free
PVC:	Free
Re use:	No
Colour:	Grey
Sterilized by:	Ethylene Oxide
Shelf life:	60 months
CFU/g	N.A.

STERILE EO



isento
Latex
free

isento
PVC
free



Nr and Revision Date:	13 -- 2019/06/27
CE Class:	Ila
CE mark:	CE0120
Directive:	MDD 93/42/EEC amended 2007/47/EC

Applicable norms:	EN ISO 13485:2016 ; EN 1041+A1:2013 ; EN ISO 11607-1:2009 ; EN ISO 11607-2:2006 ; EN ISO 11737-1:2006 ; EN ISO 11737-2:2009 ; EN 556-1:2001 ; EN ISO 14971:2012 ; EN ISO 10993-1/AC:2010 ; EN 1422:2014 ; ISO 1133-1:2011 ; EN 1041:2008 ; ENISO11737-1:2006/AC:2009 ; EN 556-1:2001/AC:2006 ; EN ISO 11135:2014 ; ISO 16142-1:2016 ; EN ISO 15 223-1:2016 ;
--------------------------	--

Level	Type	Dimension	Stamp	Labels
1	1 Rigid peel-open pack	53x173mm	stamped	
2	19 Folding box	176x056x212mm		labeled
3	76 Carton box	239 x 182 x 217		labeled

Storage instructions	Store in a dry place at room temperature and protect from sunlight. The product has a good stability and therefore, normal conditions of medical devices storage shall apply
Waste instructions	No special requirements. The waste must be disposed according to local legal regulations for the disposable of medical waste.
Intended use	Intended to be used for grasping tissues and/or holding other devices.

WASTE AND ENVIRONMENT REGULATIONS	Weight (gr) of packing material used in 100 pcs	448
--	---	-----

ISSUED BY: Quality Management Department, 2019/07/05

CONFIDENTIAL: The information contained herein is confidential and belongs to Bastos Viegas, s.a. It is provided for specific purposes and in connection with established agreements. It is not permitted disclosure, reproduction, copying or distribution, in whole or in part, for other purposes or to any third parties.