



## The optimal solution for Aspiration of Seroma

### SeromaSet – the Golden standard

- Closed system – much safer
- Easy and faster aspiration
- Less risk of contamination
- User friendly – simple to use
- Operated easily with 1 hand
- Ergonomic handle
- Patient friendly
- Latex and DEHP free
- Cost effective



### The "Old standard"

- Risk of bacterial contamination
- Not a closed system
- Bulky and time consuming
- "Hard" work for the nurse/user



### SeromaSet can be used for all Seroma aspirations, typically after:

- Breast cancer surgery – after Mastectomies, breast re-construction
- Abdominal surgery – like "Tummy Tuck" (Abdominoplasty)
- Plastic surgery – e.g. after facial surgeries
- Hernia Surgery
- C-section and Hysterectomies
- Surgery around active joints

### Ordering information:

REF	Product	Pcs per box
1060	SeromaSet, 400ml – Handle + 400ml bottle, single use	25 pcs
1061	SeromaSet, 600ml – Handle + 600ml bottle, single use	22 pcs
1020	Bottle 400ml, alone, single use	44 pcs
1021	Bottle 600ml, alone, single use	32 pcs
1050	SeromaSet Handle, alone, single use	100 pcs