




SECTION A-A




INSIDE VIEW



FRONT VIEW



SIDE VIEW



| | |
|---|------------------------|
| Device Name: Rochester-Pean Artery Forceps | |
| Size: 16cm / 160mm (6.25") | Finish: Sand |
| Code: SKG-1238-16 | Weight: 40g ±5g |
| Customer Code: | Sheet: 1 of 1 |
| Configuration: Straight & Serrated | |
| Scale: 1 : 1 | |

| | | | | |
|--|--|--------------------|-------------|--------------------|
| Material: Stainless Steel A.I.S.I 410 | Rev # | Description | Date | Approved By |
| Hardness: HRC 42 to 45 | 00 | | | |
| Dated: 10 Aug, 2017 | | | | |
| Dimensions: Metric | | | | |
| Tolerance: ± 0.50mm | Unless Specified On The Drawing Itself | | | |

| Rev # | Description | Date | Approved By |
|-------|-------------|------|-------------|
| 00 | | | |
| | | | |
| | | | |

Drawn By: aqeel-ur-rehman
CAD Designer

Approved By: Mirza Ahsan-Ul-Hussain
Q.A Manager

SURGIKARE
SURGICAL INSTRUMENTS PAKISTAN

Toor Abad, Daska Road,
Sialkot-Pakistan
URL:
www.surgikare.com