


FRONT VIEW




TIP (INSIDE)



TOTAL LENGTH 140.00MM±5.00MM

SIDE VIEW



Device Name: Thumb Tissue Forceps		Material: Stainless Steel A.I.S.I 410	Rev #	Description	Date	Approved By
Size: 14cm / 140mm (5.5")	Finish: Sand	Hardness: HRC 42 to 45	00			
Code: SKG-1140-14	Weight: 20g ±5g	Dated: 08 Aug, 2017				
Customer Code:	Sheet: 1 of 1	Dimensions: Metric				
Configuration: Straight & Serrated		Tolerance: ± 0.50mm	Unless Specified On The Drawing Itself			

Drawn By: aqeel-ur-rehman
CAD Designer

Approved By: Mirza Ahsan-Ul-Hameed
Q.A Manager

SURGIKARE
SURGERY & REPAIRS

Toor Abad, Daska Road,
Sialkot-Pakistan
URL:
www.surgikare.com