


FRONT VIEW




TIP (INSIDE VIEW)




TOTAL LENGTH 130.00MM±5.00MM

SIDE VIEW



| | | |
|--|--|---------------------|
| Device Name: Thumb Tissue Forceps | | |
| Size: 13cm / 130mm (5") | Finish: Sand | |
| Code: SKG-1140-13 | Weight: 18g ± 3g | |
| Customer Code: | Sheet: 1 of 1 | Scale: 1 : 1 |
| Configuration: Straight & Serrated | | |
| Material: Stainless Steel A.I.S.I 410 | Rev # | |
| Hardness: HRC 42 to 45 | 00 | |
| Dated: 09 Jan, 2016 | | |
| Dimensions: Metric | | |
| Tolerance: ± 0.50mm | Unless Specified On The Drawing Itself | |

Drawn By: *aqeel-ur-rehman*
CAD Designer

Approved By: 
Mr. Sami Chaudhry
Surgikare

SURGIKARE
HIGHEST QUALITY INSTRUMENTS
Toor Abad, Daska Road,
Sialkot-Pakistan
URL:
www.surgikare.com