



**FRONT VIEW**




**TIP  
INSIDE VIEW**



TOTAL LENGTH 120.00MM±5.00MM

**SIDE VIEW**



<b>Device Name:</b> Adson Tissue Forceps		<b>Material:</b> Stainless Steel A.I.S.I 410	<b>Rev #</b>	<b>Description</b>	<b>Date</b>	<b>Approved By</b>
<b>Size:</b> 12cm / 120mm (4 3/4")	<b>Finish:</b> Sand	<b>Hardness:</b> HRC 42 to 45	00			
<b>Code:</b> SKG-1153-12	<b>Weight:</b> 20g ± 5g	<b>Dated:</b> 08 Aug, 2017				
<b>Customer Code:</b>	<b>Sheet:</b> 1 of 1	<b>Scale:</b> 1 : 1				
<b>Configuration:</b> Serrated		<b>Tolerance:</b> ± 0.50mm	Unless Specified On The Drawing Itself			

**Drawn By:** *aqeel-ur-rehman*

CAD Designer

**Approved By:**



Toor Abad, Daska Road,

Sialkot-Pakistan

URL:

[www.surgikare.com](http://www.surgikare.com)