


Codman® Surgical Patties and Strips

*The benchmark  
in quality*



CODMAN Surgical Patties, sterile, X-ray detectable, 10 patties per package, 20 packages per dispenser box.

<b>80-1396</b>	5/16"	(7mm)
<b>80-1397</b>	1/4" x 1 1/2"	(7mm x 38mm)
<b>80-1398</b>	1/4" x 3"	(7mm x 76mm)
<b>80-1399</b>	1/4" x 1/4"	(7mm x 7mm)
<b>80-1400</b>	1/2" x 1/2"	(13mm x 13mm)
<b>80-1401</b>	3/4" x 3/4"	(19mm x 19mm)
<b>80-1402</b>	1/2" x 1"	(13mm x 25mm)
<b>80-1403</b>	1" x 1"	(25mm x 25mm)
<b>80-1404</b>	1/2" x 1 1/2"	(13mm x 38mm)
<b>80-1406</b>	1/2" x 2"	(13mm x 52mm)
<b>80-1407</b>	1/2" x 3"	(13mm x 76mm)
<b>80-1408</b>	1" x 3"	(25mm x 76mm)
<b>80-1409</b>	3" x 3"	(76mm x 76mm)

**CODMAN® Surgical Patties and Strips**

*CODMAN Surgical Patties and Strips provide the neurosurgeon with softness, gentleness, protection and absorbency to meet the demands of delicate intracranial and spinal surgery.*

**Absorbency**

CODMAN Patties and Strips absorb more than 5 times their weight in less than a second.

**Precision Cut**

CODMAN Patties and Strips' precision cut provides smooth edges and helps prevent linting.

**COTTONOID® Material**

COTTONOID compressed rayon provides the natural softness, absorbency and pliability of cotton. A special process interweaves individual rayon fibres to resist shedding, stretching or linting – whether wet or dry.

**Soft and Gentle**

SUTUREWELD®'s strong bond holds the indicator string to the patty so the surface which lies next to the tissue is completely smooth and free of stitches.

**Pliable and Porous**



COTTONOID material offers the advantage of being porous so CODMAN Patties and Strips easily conform to the intricate contours of delicate tissue. Yet, you can safely suction through them without aspirating the tissue or the patty strip.

**X-Ray Detectable**

Marking is where it should be – on the patty or strip.

**Unique Packaging**

A transparent, sterile envelope shows the size and shape of patties and strips, and the Count Verification Card instantly shows how many patties you started with and how many were used.



# Sponge, BICOL

**80-1460** BICOL\* Collagen Sponge, sterile.  
The BICOL Sponges are 3" (76mm) x 4" (102mm) and packaged individually sterile, sold ten packages per box.



# Strips, Neurosurgical

CODMAN Surgical Strips, sterile, X-ray detectable, 10 strips per package, 20 packages per dispenser box.

<b>80-1449</b>	1/8" x 6"	(3 x 152mm)
<b>80-1450</b>	1/4" x 6"	(6 x 152mm)
<b>80-1451</b>	1/2" x 6"	(13 x 152mm)
<b>80-1452</b>	3/4" x 6"	(19 x 152mm)
<b>80-1453</b>	1" x 6"	(25 x 152mm)
<b>80-1454</b>	1 1/2" x 6"	(38 x 152mm)
<b>80-1455</b>	2" x 6"	(51 x 152mm)
<b>80-1456</b>	3" x 6"	(76 x 152mm)
<b>80-1457</b>	3 1/2" x 6"	(89 x 152mm)



Sold by:

**MEDI PLAST**

Medioplast AB  
P.O.Box 1004

SE-212 10 Malmö, Sweden

Ph +46 40 671 23 00

mediplast.info@mediplast.com

www.mediplast.com