

M'Drain® Accessories



Ref. no. B-806805 / B-806806



Ref. no. B-806815 / B-806816



Ref. no. B-806845 / B-806846



Ref. no. B-806329



Ref. no. B-806835 / B-806836



Ref. no. B-5022021

Features:

- Kink resistant, flexible tubing
- Standard luer lock connections
- Soft, transparent tube

Connecting tubes - low pressure for CT/MR

Ref. no.	Specification	Qty/box
B-5022021	Safety patient line for CT and MRI examinations, MLL/F, 20cm. Valve to prevent backflow (Max 300 PSI)	50

Connecting tubes - low pressure

Ref. no.	Specification	Qty/box	Ref. no.	Specification	Qty/box
B-806329	Low pressure connecting tube MLL/F 20 cm (75 PSI) - sterile	25	B-806826	Low pressure connecting tube for drainage MLL/M 30 cm (150 PSI) - sterile	25
B-806805	Low pressure connecting tube for drainage bag M/Funnel 15 cm - sterile	25	B-806835	Low pressure connecting tube for drainage 15 cm with 3 ways stopcock M/Funnel - sterile	25
B-806806	Low pressure connecting tube for drainage bag M/Funnel 30 cm - sterile	25	B-806836	Low pressure connecting tube for drainage 30 cm with 3 ways stopcock M/Funnel - sterile	25
B-806815	Low pressure connecting tube for drainage bag M/Funnel 15 cm with 1 way stopcock - sterile	25	B-806845	Low pressure connecting tube for drainage MLL/15 cm W/3W. Mills off stopcock (75 PSI) - sterile	25
B-806816	Low pressure connecting tube for drainage bag M/Funnel 30 cm with 1 way stopcock - sterile	25	B-806846	Low pressure connecting tube for drainage MLL/30 cm W/3W. Mills off stopcock (75 PSI) - sterile	25

DS-2-01 Version 10

Contact/Place order: Mermaid Medical A/S · Frydensbergvej 25 · DK-3660 Stenløse · www.mermaidmedical.com · customer@mermaidmedical.com



Denmark

P +45 47 10 85 70
F +45 47 10 85 71

Sweden

P +46 31 755 41 28
F +46 31 755 45 99

Finland

P +358 923 113238
F +358 923 113239

Norway/Iceland

P +47 21 53 88 94
F +47 21 53 88 95

United Kingdom

P +44(0)1704 560493
F +44(0)1704 551182

Spain

P +34 9253 33159
F +34 9253 33159
customer@mermaidmedical.es


Netherlands

P +31 314 63 2070
F +31 314 63 1685
info@mermaidmedical.nl

Germany

P +49 40 8518665-0
F +49 40 8518665-1
info@mermaidmedical.de

M'Drain® Accessories



A10500750



A10602000C



A10500350

Features:

- Flocked back reduces skin discomfort
- Drainage bag with anti-slip leg straps
- Anti-reflux valve prevents retrograde bacterial contamination
- Discrete size - invisible beneath clothes
- Anti-reflux valve minimizes risk of infection

Drainage bag

Ref. no.	Specification	Qty/box
A10500350	Drainage bag - 350 ml - female - sterile incl. 5 pcs. anti-slip leg straps	50
A10500750	Drainage bag - 750 ml - female - sterile incl. 2 pcs. anti-slip leg straps	50
A10602000C	Drainage bag - 2000 ml - male - sterile incl. 5 pcs. anti-slip leg straps with Clamp	50

DS-2-01 Version 10



Contact/Place order: Mermaid Medical A/S · Frydensbergvej 25 · DK-3660 Stenløse · www.mermaidmedical.com · customer@mermaidmedical.com

Denmark
P +45 47 10 85 70
F +45 47 10 85 71

Sweden
P +46 31 755 41 28
F +46 31 755 45 99

Finland
P +358 923 113238
F +358 923 113239

Norway/Iceland
P +47 21 53 88 94
F +47 21 53 88 95

United Kingdom
P +44(0)1704 560493
F +44(0)1704 551182

Spain
P +34 9253 33159
F +34 9253 33159
customer@mermaidmedical.es

Netherlands
P +31 314 63 2070
F +31 314 63 1685
info@mermaidmedical.nl

Germany
P +49 40 8518665-0
F +49 40 8518665-1
info@mermaidmedical.de