



Product description

| | |
|---------------------|--|
| Order number | 95.963.005 |
| Product description | Mounting plate (max. width: 10 mm) for system bed rails for combination with the wall holder euroMatic® (95.963.001) |

Product characteristics

| | |
|--------------|--|
| Design | Mounting plate for the wall holder euroMatic® (95.963.001) |
| Application | For attachment to system bed rails in combination with the wall holder euroMatic® (95.963.001) |
| Suitable for | Wall holder euroMatic® (95.963.001) |

Size

| | |
|-------------------|--------|
| Width of product | 15 mm |
| Height of product | 79 mm |
| Length of product | 130 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Stainless steel

Purity & certification

Product in accordance with MDD or IVDD No marking necessary

Lot or serial number non batched

Packaging

| | |
|--------------------|---------------|
| Minimum order qty. | 1 |
| Pcs. / inner box | 1 |
| Pcs. / box | 20 |
| Pcs. / pallet | 6480 |
| Depth of case | 268 mm |
| Width of case | 191 mm |
| Height of case | 100 mm |
| Case volume | 0,005 cbm |
| Weight of case | 3,68 kg |
| EAN of inner box | 4038917364096 |
| EAN of outer box | 4038917364089 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

