



TVÄTTSET

ART.NR 60171116T

	Beskrivning	Avd/Transport Box
1	Skål 1000 ml, blå	6/36
8	Tork, Nonwoven Ø 50 mm	
1	Peang blå, 12 cm	

Material:

Bowls: Polypropylene PP

Peang: Polypropylene PP

Balls/Swabs: 70 % Rayon 30 % Polyester Peroxide bleached

Packaging: Peel-pack, Paper/PE

Sterilization: Shelf-life 5 years from production date

Manufacturer: Mediplast AB

MDD classification: Is

To be stored in room temperature

Transport Box = minimum sales unit

Destruction: Combustion, recyclable carton, recyclable well-cardboard



STERILE EO



Mediplast AB

P.O. Box 1004

SE-212 10 Malmö, Sweden

Bronsåldersgatan 2, SE-213 76 Malmö, Sweden

Ph. no. +46 (0)40 671 23 00

Fax no. +46 (0)40 671 23 10

Customer service (Sweden only) 020 788035

www.mediplast.com
Org. nr. 556140-9011