



## TVÄTTSET

### ART.NR 60171022T

	Beskrivning	Avd/Transport Box
1	Omslagsduk 70 x 70 cm	4/16
1	Rondskål 700 ml, blå	
1	Skål, 500 ml, blå	
8	Tork, Nonwoven Ø 50 mm	
1	Peang m ring, 18 cm	

#### Material:

Fields: 100 % Cellulose tissue + PE Polyethylene  
 Bowls: Polypropylene PP  
 Peang: Acrylonitrile Butadiene Styrene ABS  
 Balls/Swabs: 70 % Rayon 30 % Polyester Peroxide bleached  
 Packaging: Peel-pack, Paper/PE  
 Sterilization: Shelf-life 5 years from production date  
 Manufacturer: Mediplast AB  
 MDD classification: Is  
 To be stored in room temperature  
 Transport Box = minimum sales unit  
 Destruction: Combustion, recyclable carton, recyclable well-cardboard



**STERILE EO**



#### Mediplast AB

P.O. Box 1004  
 SE-212 10 Malmö, Sweden  
 Bronsåldersgatan 2, SE-213 76 Malmö, Sweden  
 Ph. no. +46 (0)40 671 23 00  
 Fax no. +46 (0)40 671 23 10  
 Customer service (Sweden only) 020 788035

[www.mediplast.com](http://www.mediplast.com)  
 Org. nr. 556140-9011