

Water sampling bottles, sterile



Leverantör: VWR Collection

PET, virgin, food grade with PP caps.



- ▶ Different coloured caps for easy identification of test
- ▶ Tamper-evident strip to ensure sampling integrity
- ▶ Full traceability, batch number and individual bottle number on each label
- ▶ Fully recyclable
- ▶ Cap diameter: 32 mm

Lightweight, durable bottles. Caps designed to shield neck from contamination. Choice of dosed and undosed bottles. Rectangular bottles for better storage with limited space.







Blue cap: For natural, non chlorinated water without sodium thiosulphate

Red cap: For chlorinated water dosed with 20 mg/l sodium thiosulphate

White cap: For chlorinated water dosed with 120 mg/l sodium thiosulphate

ORDER

Square bottle

VWR art.nr	Antal	Volym	Dos	BxDxH	Lock, färg	Förp.
331-0268		500 ml	Ingen	70x70x192 mm	Blå	108
331-0272		500 ml	Natriumtiosulfat (20 mg/liter)	70x70x192 mm	Röd	108
331-0276		500 ml	Natriumtiosulfat (120 mg/liter)	70x70x192 mm	Vit	108
331-0269		1000 ml	Ingen	89x89x222 mm	Blå	70
331-0273		1000 ml	Natriumtiosulfat (20 mg/liter)	89x89x222 mm	Röd	70
331-0277		1000 ml	Natriumtiosulfat (120 mg/liter)	89x89x222 mm	Vit	70

Rectangular Bottle

VWR art.nr	Antal	Volym	Dos	BxDxH	Lock, färg	Förp.
------------	-------	-------	-----	-------	------------	-------

331-0256		350 ml	Ingen	70x50x174 mm	Blå	144
331-0270		350 ml	Natriumtiosulfat (20 mg/liter)	70x50x174 mm	Röd	144
331-0274		350 ml	Natriumtiosulfat (120 mg/liter)	70x50x174 mm	Vit	144
331-0267		500 ml	Ingen	70x50x212 mm	Blå	144
331-0271		500 ml	Natriumtiosulfat (20 mg/liter)	70x50x212 mm	Röd	144
331-0275		500 ml	Natriumtiosulfat (120 mg/liter)	70x50x212 mm	Vit	144