



## TVÄTTSET ART.NR 605158T

	Beskrivning	Avd/Transport Box
1	Omslagsduk 50 x 70 cm	8/32
1	Rondskål 700 ml, vit	
1	Skål 500 ml	
8	Kompresser 7,5 x 7,5 cm Nonwoven	
8	Tork Ø 50 mm Nonwoven	
1	Peang m ring, 18 cm	

**Material:**

Wrapping field: 100% Cellulose tissue + Polyethylene PE

Kidney Dish: Polystyrene PS

Bowl: Polypropylene PP

Balls/Swabs: Paper 70% Rayon 30 % Polyester

Peang: Acrylonitrile Butadiene Styrene ABS

Packaging: Peel-pack Paper/PE/PA

Sterilization: Shelf-life 5 years from production date

Manufacturer: Medioplast AB

MDD Classification: Is

To be stored in room temperature

Transport Box = minimum sales unit

Destruction: Combustion, recyclable carton, recyclable well-cardboard



**STERILE EO**



**Medioplast AB**

P.O. Box 9504

SE-200 39 Malmö, Sweden

Kantyxegatan 29, SE-213 76 Malmö, Sweden

Ph. no. +46 (0)40 671 23 00

Fax no. +46 (0)40 671 23 10

Customer service (Sweden only) 020 788035

www.mediplast.com  
Org. nr. 556140-9011