

Ordering information

1000 μ l PIPETTE TIPS



PROVEN &
DIVERSE

Universal

The versatile 1,000 µl tip can be used with the Bag, Box, SingleRefill or StackPack. This is also available in the exceptionally high purity Biosphere® plus and in the low retention version!

Reliable **DE-M**

The proven tip design in our new, innovative box system combines well-known design and revolutionary handling. The precision and reliability of our tips conform to the international test standard ISO 8655-2, making them the perfect partner for your pipette.

Extended range

The 1,000 µl tip is now also available in our StackPack design. This allows for fast processing directly from the StackPack and uncomplicated refilling. Is special protection and the highest purity of tips of high importance to you? Then try out our SingleRefill.



TOP TIP Identification at a glance! In addition to the familiar colour code of the trays, you will now also find the most important information about the tip directly on each tray.

Ordering information

Vol.	Order no.	Description	Filter	Purity	Design	Packaging	Pieces/ carton
1,000 µl	70.3050.275	Low retention filter tip, 1000 µl, Biosphere® plus, transparent, 96 pcs/Box	+		Low Retention		1,920
	70.3050.375	Low retention filter tip, 1000 µl, Biosphere® plus, transparent, 96 pcs/SingleRefill	+		Low Retention		1,920
	70.3050.255	Filter tip, 1000 µl, Biosphere® plus, transparent, 96 pcs/Box	+		Standard		1,920
	70.3050.355	Filter tip, 1000 µl, Biosphere® plus, transparent, 96 pcs/SingleRefill	+		Standard		1,920
	70.3050.205	Pipette tip, 1000 µl, Biosphere® plus, transparent, 96 pcs/Box	-		Standard		1,920
	70.3050.305	Pipette tip, 1000 µl, Biosphere® plus, transparent, 96 pcs/SingleRefill	-		Standard		1,920
	70.3050.210	Low retention pipette tip, 1000 µl, PCR Performance Tested, transparent, 96 pcs/Box	-		Low Retention		1,920
	70.3050.110	Low retention pipette tip, 1000 µl, PCR Performance Tested, transparent, 5x 96 pcs/StackPack	-		Low Retention		1,920
	70.3050.200	Pipette tip, 1000 µl, PCR Performance Tested, transparent, 96 pcs/Box	-		Standard		1,920
	70.3050.100	Pipette tip, 1000 µl, PCR Performance Tested, transparent, 5x 96 pcs/StackPack	-		Standard		1,920
	70.3050	Pipette tip, 1000 µl, PCR Performance Tested, transparent, 500 pcs/Bag	-		Standard		5,000
	70.3050.020	Pipette tip, 1000 µl, PCR Performance Tested, blue, 500 pcs/Bag	-		Colour variation		5,000
	70.3050.030	Conductive pipette tip, 1000 µl, PCR Performance Tested, black, 500 pcs/Bag	-		Conductive		5,000

We print on FSC certified paper.



Further information can be found at www.refillrevolution.tips

SARSTEDT AG & Co. KG
P.O. Box 12 20
D-51582 Nümbrecht
Phone: +49 2293 305 0
Fax: +49 2293 305 3992
export@sarstedt.com
www.sarstedt.com