

Ordering information

300 μ l PIPETTE TIPS



LITHE & AGILE

Ergonomic

Our new extra slim tip can enter any reaction vessel, screw cap micro tube and small tube with ease without even touching the vessel wall. Thanks to its 60-mm length, it is able to reach the bottom of the vessel making safe pipetting even easier. The volume level can be checked with the help of the clear level rings.

Innovative

The sophisticated, new packaging system means that the tips can be used for many different applications. Our system adapts to your needs. Whether high-purity work with the SingleRefill, high throughput with the StackPack or simply a filled Box – there is a suitable solution for every requirement.

Perfect fit

Our tips fit perfectly on all standard pipettes. Have a look at www.sarstedt.com to see which tip fits your pipette!



TOP TIP Do you work with particularly viscous or expensive substances?
 ■ Then test our low retention tips which provide even better sample recovery!

Ordering information

Vol.	Order no.	Description	Filter	Purity	Design	Packaging	Pieces/carton
300 µl	70.3040.275	Low retention filter tip, 300 µl, Biosphere® plus, transparent, 96 pcs/Box	+		Low retention		1,920
	70.3040.375	Low retention filter tip, 300 µl, Biosphere® plus, transparent, 96 pcs/SingleRefill	+		Low retention		1,920
	70.3040.255	Filter tip, 300 µl, Biosphere® plus, transparent, 96 pcs/Box	+		Standard		1,920
	70.3040.355	Filter tip, 300 µl, Biosphere® plus, transparent, 96 pcs/SingleRefill	+		Standard		1,920
	70.3040.205	Pipette tip, 300 µl, Biosphere® plus, transparent, 96 pcs/Box	-		Standard		1,920
	70.3040.305	Pipette tip, 300 µl, Biosphere® plus, transparent, 96 pcs/SingleRefill	-		Standard		1,920
	70.3040.210	Low retention pipette tip, 300 µl, PCR Performance Tested, transparent, 96 pcs/Box	-		Low retention		1,920
	70.3040.110	Low retention pipette tip, 300 µl, PCR Performance Tested, transparent, 5x 96 pcs/StackPack	-		Low retention		1,920
	70.3040.200	Pipette tip, 300 µl, PCR Performance Tested, transparent, 96 pcs/Box	-		Standard		1,920
	70.3040.100	Pipette tip, 300 µl, PCR Performance Tested, transparent, 5x 96 pcs/StackPack	-		Standard		1,920
	70.3040	Pipette tip, 300 µl, PCR Performance Tested, transparent, 1000 pcs/Bag	-		Standard		10,000

We print on FSC certified paper.



Further information can be found at www.refillrevolution.tips

SARSTEDT AG & Co. KG
 P.O. Box 12 20
 D-51582 Nümbrecht
 Phone: +49 2293 305 0
 Fax: +49 2293 305 3992
export@sarstedt.com
www.sarstedt.com

SARSTEDT