

ETHILON*

(Polyamide 6) Suture



ETHILON*

Monofilament Polyamide 6

SUTURE

Description: **Monofilament** For smooth passage through tissue.
 Synthetic For minimal tissue reaction.
 Non-absorbable Provides prolonged tensile strength retention in tissue.

Color: Blue, black or clear.

Range: 11-0 to 2 (USP). Supplied as needled sutures armed with, for example, MULTIPASS* and PRIME* needles. ETHILON* is also available as ligatures.

Indications: ETHILON* (black and clear) sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, microsurgical and neurosurgical procedures. Blue ETHILON* sutures are intended for use in skin closure. Typical areas of use include:











- Cuticular sutures
- Nerve adaptation
- Ophthalmology

Full details in the Instruction for Use included in every package.












ETHILON* Monofilament Polyamide 6

3/8 Circle Reverse Cutting Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 C-2	13 mm	6-0	BLACK	45	250 microns	36	667H
		5-0	BLACK	45	250 microns	36	668H
 M-2	13 mm	5-0	BLACK	45	360 microns	36	644H
 FS-3	16 mm	6-0	BLACK	45		36	660H
		6-0	BLACK	45	VISI-BLACK*	36	EH7177H
		5-0	BLACK	45		36	1628H
		4-0	BLACK	45		36	1162H
 FS-2	19 mm	5-0	BLACK	45		36	661H
		4-0	BLUE	45		36	EH7144H
		4-0	BLACK	75		36	EH7143BH
		3-0	BLACK	45		36	653H
		2-0	BLUE	90		36	F3217
 FS-2S	19 mm	4-0	BLACK	45	Slimmed	36	662SLH
 FS-1	24 mm	5-0	BLACK	45		36	EH7640H
		4-0	BLACK	45		36	1629H
		4-0	BLACK	75		36	FJ5272BH
		3-0	BLACK	45		36	663H
		3-0	BLACK	75		36	669H
		3-0	BLUE	90		36	F3221
		2-0	BLACK	45		36	EH7826BH
		2-0	BLUE	75		36	EH7827H
 FS	26 mm	4-0	BLACK	45		12	W345BG
		3-0	BLACK	45		12	W320BG
		2-0	BLACK	45		36	664H
 FSL	30 mm	4-0	BLACK	75		36	1670BH
		3-0	BLACK	45		36	1171H
		3-0	BLACK	75		36	1671H
		3-0	BLACK	45		36	S2455
		2-0	BLUE	75		36	EH7798H
 FSLX	36 mm	3-0	BLACK	75		36	1673BH
		2-0	BLUE	75		36	1674H
		0	BLACK	75		36	P7085H
 FSLY	45 mm	2-0	BLACK	100		12	W736G






ETHILON* Monofilament Polyamide 6

3/8 Circle Reverse Cutting PRIME* Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 P-6	8 mm	7-0	BLACK	45		36	1647H
 P-1	11 mm	7-0	BLACK	45		36	1696H
		7-0	BLACK	45	MULTIPASS*	36	MPE1696H
		6-0	BLACK	45		36	697H
		6-0	BLACK	45	MULTIPASS*	36	MPE697H
		5-0	BLACK	45		24	W1611T
 P-3	13 mm	6-0	BLACK	45	MULTIPASS*	36	1698H
		5-0	BLACK	45	MULTIPASS*	36	698H
		4-0	UNDYED	45		36	691H
		4-0	BLACK	45	MULTIPASS*	36	699H
 PS-3	16 mm	6-0	BLACK	45	MULTIPASS*	12	1665G
		5-0	BLACK	45	MULTIPASS*	36	1668H
		4-0	BLACK	45	MULTIPASS*	36	EH7761H
 PS-2	19 mm	5-0	BLACK	45	MULTIPASS*	36	1666H
		4-0	BLUE	45		36	EH7162H
		4-0	BLUE	45		24	W1620T
		4-0	BLUE	45	MULTIPASS*	36	MPE1855H
		4-0	UNDYED	45		12	1611G
		3-0	BLACK	45	MULTIPASS*	36	1669H
 PS-2 S	19 mm	4-0	BLACK	45	MULTIPASS*	36	1667SLH
 PS-1	24 mm	4-0	BLACK	75	MULTIPASS*	36	EH7957BH
		3-0	BLACK	45		36	1663H
		3-0	BLACK	75	MULTIPASS*	36	EH7563H
 PS	26 mm	3-0	BLUE	45	MULTIPASS*	36	MPE1625H
		3-0	BLACK	45	MULTIPASS*	36	W1625BH
		3-0	BLACK	75	MULTIPASS*	36	W1685BH
		2-0	BLUE	75		24	W1626T
 PSL	30 mm	3-0	BLUE	75		36	EH7113H

ETHILON* Monofilament Polyamide 6


3/8 Circle Conventional Cutting PRIME* Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
	13 mm	6-0	BLACK	45		36	JM5140H
	16 mm	6-0	BLACK	45	MULTIPASS*	12	1866G
		5-0	BLACK	45	MULTIPASS*	36	1865H
	19 mm	4-0	BLACK	45	MULTIPASS*	12	1894G
		1	BLACK	2x50	With retention tubing	12	W2797
		1	BLACK	75		12	W689

1/2 Circle Reverse Cutting Needles



Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
	13 mm	5-0	BLACK	45		36	648H
	17 mm	5-0	BLACK	75		36	F2423H
	22 mm	3-0	BLACK	45		12	642G
	48 mm	1	BLACK	100		12	W738

Straight Cutting Needles














Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
	3-0	BLACK	75		36	627H
	2-0	BLACK	75		36	628H

ETHILON* Monofilament Polyamide 6

1/2 Circle Taper Point Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
	40 mm	0	BLACK	150	Loop	12	W740
	48 mm	1	BLACK	150	Loop	12	W748

3/8 Circle Taper Point Needles




Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
BV130-3 	3,7 mm, 150 microns	9-0	BLACK	13		12	FG2819
BV75-3 	3,8 mm, 75 microns	10-0	BLACK	13		12	W2870
		9-0	BLACK	13		12	W2871
BV50-3 	3,8 mm, 50 microns	11-0	BLACK	13		12	W2881
BV130-4 	4,7 mm, 150 microns	9-0	BLACK	13		12	W2813
		8-0	BLACK	13		12	W2812
BVH100-3 	5 mm, 100 microns	10-0	BLACK	13	1/2 circle	12	W2860
BV75-4 	5,1 mm, 75 microns	10-0	BLACK	13		12	W2850
BV-4 	5,1 mm, 100 microns	10-0	BLACK	15		12	2814G
BV100-4 	5,1 mm, 100 microns	10-0	BLACK	13		12	W2830
		9-0	BLACK	13		12	W2829
		9-0	BLACK	15		12	EH7448G
BV-2 	6,5 mm, 150 microns	10-0	BLACK	15		12	2810G
		9-0	BLACK	15		12	2809G
		8-0	BLACK	15		12	2808G
BV130-5 	6,5 mm, 130 microns	10-0	BLACK	13		12	W2810
		8-0	BLACK	13		12	W2808
BV130-5 	6,5 mm, 130 microns	10-0	BLACK	30		12	W2801
BV-1  BV-1	9,3 mm, 250 microns	8-0	BLACK	45		12	W8170
BV-1 	9,3 mm, 250 microns	7-0	BLACK	45		36	EH7446H

ETHILON* Monofilament Polyamide 6


3/8 Circle CS ULTIMA* Spatula

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-140-6 	6,5 mm, 150 microns	10-0	BLACK	15		12	7060G
		10-0	BLACK	30		12	U7061
CS-140-6 	6,5 mm, 150 microns	10-0	BLACK	30		12	U7003



3/8 Circle MICROPOINT* Spatula

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-35-C 	6 mm, 150 microns	10-0	BLACK	15		12	EH7980G
GS-9 	6,5 mm, 200 microns	10-0	BLACK	30		12	1718G
		9-0	BLACK	30		12	1717G
TG140-6  TG140-6	6,5 mm, 150 microns	11-0	BLACK	30		12	W1780
		10-0	BLACK	30		12	W1770
		9-0	BLACK	30		12	W1769


1/2 Circle MICROPOINT* Reverse Cutting Needles


Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
G-7 	8 mm	8-0	BLACK	45		12	W1765



1/2 Circle MICROPOINT* Spatula



Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
TG160-6 	5,5 mm, 150 microns	10-0	BLACK	30		12	W1756
TG175-8 	7 mm, 200 microns	8-0	BLACK	30		12	1716G

1/2 Circle CS ULTIMA* Spatula

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-160-6 	5,5 mm, 150 microns	10-0	BLACK	30		12	U7000

1/4 Circle MICROPOINT* Spatula							
		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
TG100-8 	6,5 mm, 220 microns	8-0	BLACK	30		12	B7714

ETHILON* Monofilament Polyamide 6							
1/4 Circle Spatula Needles							
		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
S-14  S-14 	8 mm	5-0	BLACK	30		12	B7712

1/4 Circle TAPERCUT* Needles							
		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
V100-3 	3,5 mm, 100 microns	10-0	BLACK	13		12	W2899
		9-0	BLACK	13		12	W2898
V150-4 	4,5 mm, 150 microns	10-0	BLACK	15		12	W1701

Ligatures							
		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
LIGAPAK*		2-0	BLACK	250	Reel	20	F4354