


DISPOSABLE BIOPSY PUNCH WITH PLUNGER SYSTEM



STERILE • STAINLESS STEEL • SINGLE USE • COLOR CODED • PATENT PENDING



Seamless Stainless Steel Cutting Edge
Ensures a consistently smooth cutting surface, allowing optimal specimen quality



Internal Plunger System

Ejects lodged skin specimen inside the metal lumen of the punch



NEW	BPP-10F 1.0mm w/plunger	
NEW	BPP-15F 1.5mm w/plunger	
NEW	BPP-20F 2.0mm w/plunger	
NEW	BPP-30F 3.0mm w/plunger	
NEW	BPP-40F 4.0mm w/plunger	

Punches are available in three sizes, marked clearly on the handle, blister pack and box
Punches are individually packed in color-coded, transparent, sturdy blister packs, 20 per dispenser box

Product specifications may be changed without notice



CE
0197



kai industries co.,ltd.
1110, Oyana, Seki City, Gifu Pref. 501-3992 Japan
Phone : +81 (0) 575 28 8600 Fax : +81 (0) 575 28 8611
Email: info-ind@kai-group.com http://www.kai-ind.co.jp



kai Europe GmbH
Kottendorfer Str. 5, 42967 Solingen, Germany
Phone : +49 (0) 212 232380 Fax : +49 (0) 212 232389

KIM-020 0510/ 0611 0,1K