

Pulmodyne®

BiTrac Select MaxShield™

ONE Mask. COMPLETE Coverage.

Pulmodyne's BiTrac Select MaxShield™ is a full coverage NIV mask with multi-purpose interchangeable elbows. The full facial perimeter can simplify mask fitting, reduce nasal skin breakdown, and increase patient comfort.



Standard



Anti-Asphyxia
with Leak



Standard
Bronchoscopy



For more information about BiTrac Select MaxShield™ and our complete line of ventilation products visit our website at www.pulmodyne.com or call 317.246.5505.

....bringing change to life®



BiTrac MaxShield Select™



BiTrac SE MaxShield Select™



BiTrac™ NIV Circuit

Part No.	Size	Case/Qty.
----------	------	-----------

Masks (All masks include Head Strap)

BiTrac MaxShield Select™ with Interchangeable Anti-Asphyxia Elbow with Leak 22mm Female, and MaxShield Head Strap

313-9555	Extra Large	10 per case
313-9552	Adult	10 per case
313-9554	Small Adult	10 per case

BiTrac SE MaxShield Select™ with Interchangeable Standard (SE) Elbow 22mm Female, and MaxShield Head Strap

313-9556	Extra Large	10 per case
313-9553	Adult	10 per case
313-9551	Small Adult	10 per case

Interchangeable Elbow Options

313-9599	Standard (SE) Elbow - for dual limb circuit	10 per case
313-9593	Anti-Asphyxia Elbow with Leak	10 per case
313-9597	Standard (SE) Bronchoscopy Elbow - for dual limb circuit	10 per case

Circuits

BiTrac™ NIV Circuit (Bilevel/CPAP Applications) 72" Smooth ID Tubing w/ swivel connector w/ exhalation port, cap plug; includes optional 24" tubing w/ "T" and proximal pressure line

313-9074	10 per case
----------	-------------

BiTrac™ NIV Circuit (Bilevel/CPAP Applications) 72" Smooth ID Tubing w/ attached filter, swivel connector w/ exhalation port, cap plug; includes optional 24" tubing w/ "T" and proximal pressure line

313-9075	10 per case
----------	-------------

ValueTrac™ NIV Circuit (Bilevel/CPAP Applications) 72" Corrugated Tubing and swivel connector w/ exhalation port and proximal pressure line

313-9078	10 per case
----------	-------------



Made in USA

Single Patient Use

Rx ONLY



Latex Free

Pulmonary
2055 Executive Drive
Indianapolis, IN 46241 USA

T: 317.246.5505
F: 317.246.5501

www.pulmonary.com

Pulmonary®

bringing change to life®

"Pulmonary" and "bringing change to life" are registered trademarks of Pulmonary, Inc. BiTrac, BiTrac MaxShield, BiTrac SE MaxShield, and ValueTrac are trademarks of Pulmonary, Inc.