

SANDEL[®] ERGO-STEP[™] STOOL

STAFF & PATIENT ERGONOMICS

Light-weight, stackable step stools

SANDEL[®] Ergo-Step[™] stools are light-weight and stackable to provide a safe, comfortable working platform. Ergo-Step stools feature a cushioned top to help prevent foot and leg fatigue during long procedures.

FEATURES AND BENEFITS:

- Lightweight - Approximately 7 lbs.
- Anti-skid feet keep stools firmly in place
- Offered in top and bottom versions that can be safely stacked
- Interlocking feature means stools can be attached side-by-side to create a safe working platform
- Weight capacity: 500 lbs (227 kg.)
- For added comfort, a 13" x 17" SANDEL[®] ErgoPlus[™] Anti-Fatigue Mat fits neatly on top
- Ergo-Step[™] Cart available for easy transport of stools
- Bright orange and highly recognizable

NOTE: Top step stools (1170-T) should only be stacked on bottom step stools (1170-B). 1170-B can be stacked interchangeably as either a top or a bottom stool.

