


<b>Product Name</b>	Alterna™, ClearFlex™ economy silicone, anaesthetic face mask, size 2, paediatric, no hook ring, 22F	
<b>Reference numbers</b>	8848002	
<b>Product Type</b>	Alterna Anaesthetic Mask	
<b>Manufacturer</b>	Intersurgical	
<b>Materials</b>	 Si	
<b>Packaging</b>		
<b>Unit Packaging</b>	Individually packaged	
<b>Box quantity</b>	1	
<b>Pallet quantity</b>	320	
<b>Instructions for Use</b>	Instructions for use with product	Yes
	Instructions for use with each box	No
<b>CE Class</b>	IIa	
<b>Applicable Standards</b>	BS EN ISO 13485: Medical devices — Quality management systems — Requirements for regulatory purposes ISO 9001: Quality Management Systems - Requirements ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied BS ISO 18190: Anaesthetic and respiratory equipment - General requirements for airways and related equipment	
<b>Sterile</b>	No	
<b>Disposal Methods</b>	As hospital guidelines	
<b>Re-use</b>	Reprocessible - 50 cycles	

## Product characteristics

### Technical information

Patient size	Paediatric
Size	2
Hook ring	No
Approximate fill volume (ml)	78.0
Weight (g)	39.1
Connection	22F elastomeric



**Material  
Abbreviations**

ABS	Acrylonitrile Butadiene Styrene
AL	Aluminium
CR	Polychloroprene Rubber
EDPM	Ethylene Propylene Terpolymer Rubber
HDPE	High Density Polyethylene
HIPS	High Impact Polystyrene
LDPE	Low Density Polyethylene
LSR	Liquid Silicone Rubber
MBS	Methacrylate Butadiene Styrene
NRL	Natural Rubber Latex
PA	Polyamide (Nylon)
PC	Polycarbonate
PE	Polyethylene
POM	Polyoxymethylene
PP	Polypropylene
PS	Polystyrene
PSU	Polysulfone
PU	Polyurethane
PVC	Poly Vinyl Chloride
SB	Styrene Butadiene
Si	Silicone
SS	Stainless Steel
TPE	Thermoplastic Elastomer
SBC	Styrene Butadiene Copolymer
SEBS	Styrene-Ethylene-Butadiene-Styrene

