



Product description

Order number	70.3031.255
Product description	Filter tip, working volume: 200 µl, transparent, filling level rings, Biosphere® plus, suitable for Eppendorf, Gilson, Finnpiquette, Biohit, Brand and Socorex and products of identical design, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpiquette, Biohit, Brand and Socorex and products of identical design
Volume of work	200 µl
Filter	yes
Coding colour	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Max. volume	200 µl
Length of product	51 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Filter material	Polyethylene (PE)

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	1920
Type of inner packaging	box
Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	96000
Depth of case	348 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0,0277 cbm
Weight of product	0,0011 kg
Weight of case	5,16 kg
EAN of inner box	4038917794091
EAN case	4038917794107

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

