

Functional Human Knee Joint Model with Ligaments - 3B Smart Anatomy

Product Information



| | |
|------------|-----------------|
| Weight | 0.705 kg |
| Dimensions | 12 x 12 x 34 cm |
| Brand | 3B Scientific |
| Item No. | A82 [1000163] |

- This product carries a 5 year warranty

Description

This high quality, life-size functional knee [joint model](#) clearly shows the anatomy and mechanics of the [knee joint](#). This fully flexible knee joint model demonstrates abduction, anteversion, retroversion and internal/external rotation.

Model consists of portion of [femur, tibia and portion of fibula](#); also includes meniscus, patella with quadriceps tendon and joint ligaments, including the ACL and PCL.

Delivered on removable stand for easy study or display.

https://www.3bscientific.com/se/functional-human-knee-joint-model-with-ligaments-3b-smart-anatomy.p_35_160.html

