



Personalized PDF Catalog
Catalogue Generated den 26 april 2022

Pipette tips, epT.I.P.S.® reloads (General Lab Product)

Leverantör: EPPENDORF



A compact system for refilling epT.I.P.S.® Available as one filled box of tips; a set comprising one box and five trays of tips; and reloads of ten tip trays. Components are 100% recyclable.

- Optimised system for multi-channel pipettes
- Reloads available in a PCR clean version: Free of human DNA, DNase, RNase and PCR inhibitors

epT.I.P.S., in combination with adequate pipettes, are intended for dispensing liquids in the volume range of 0,1 µl to 10 ml.

In vivo applications (applications in or on the human body) are not allowed. epT.I.P.S. and the appropriate pipettes only must be used by correspondingly trained specialists.

The user must have read thoroughly the operating instructions of the pipette and the instructions for use of the epT.I.P.S. and must have made himself familiar with the functioning of the pipette with the pipette tips.


ORDER

epT.I.P.S. reload



VWR art.nr	Antal	Volym	Längd	Färgkod	Förp.	Förp.
613-6068		0,1 - 10 µl	34 mm	Mörkgrå	10 trays, 96 each	960
613-6069		0,1 - 20 µl	40 mm	Grå	10 trays, 96 each	960
EPPE0030076.044		0,1 - 20 µl	42 mm	Vit	10 trays, 384 each	3.840
613-6070		0,5 - 20 µl	46 mm	Ljusgrå	10 trays, 96 each	960
EPPE0030076.052		0,5 - 100 µl	53 mm	Ljusgul	10 trays, 384 each	3.840
613-6071		2 - 200 µl	53 mm	Gul	10 trays, 96 each	960
613-6072		20 - 300 µl	55 mm	Orange	10 trays, 96 each	960

613-6073		50 - 1000 µl	71 mm	Blå	10 trays, 96 each	960
613-6074		50 - 1250 µl	76 mm	Grön	10 trays, 96 each	960
613-6076		50 - 1250 µl	103 mm	Mörkgrön	10 trays, 96 each	960
613-6075		250 - 2500 µl	115 mm	Röd	10 trays, 48 each	480

epT.I.P.S. reloads PCR clea

VWR art.nr	Antal	Volym	Längd	Färgkod	Förp.	Förp.
EPPE0030073.754		0,1 - 10 µl	34 mm	Mörkgrå	10 trays, 96 each	960
EPPE0030073.770		0,1 - 20 µl	40 mm	Grå	10 trays, 96 each	960
EPPE0030076.001		0,1 - 20 µl	42 mm	Vit	10 trays, 384 each	3.840
EPPE0030073.797		0,5 - 20 µl	46 mm	Ljusgrå	10 trays, 96 each	960
EPPE0030076.010		0,5 - 100 µl	53 mm	Ljusbgul	10 trays, 384 each	3.840
EPPE0030073.819		2 - 200 µl	53 mm	Gul	10 trays, 96 each	960
EPPE0030073.835		20 - 300 µl	55 mm	Orange	10 trays, 96 each	960
EPPE0030073.851		50 - 1000 µl	71 mm	Blå	10 trays, 96 each	960
EPPE0030073.878		50 - 1250 µl	76 mm	Grön	10 trays, 96 each	960
EPPE0030073.649		50 - 1250 µl	103 mm	Mörkgrön	10 trays, 96 each	960
EPPE0030073.894		250 - 2500 µl	115 mm	Röd	10 trays, 48 each	480

epT.I.P.S. bo

VWR art.nr	Antal	Volym	Längd	Färgkod	Förp.	Förp.
EPPE0030073.010		0,1 - 10 µl	34 mm	Mörkgrå	1 box incl. 96 tips	96
EPPE0030073.150		250 - 2500 µl	115 mm	Röd	1 box incl. 48 tips	48

epT.I.P.S. se

VWR art.nr	Antal	Volym	Längd	Färgkod	Förp.	Förp.
------------	-------	-------	-------	---------	-------	-------







EPPE0030073.215		0,1 - 10 µl	34 mm	Mörkgrå	1 box incl. 5 trays, 96 each	480
EPPE0030076.028		0,1 - 20 µl	42 mm	Vit	1 box incl. 5 trays, 384 each	1.920
EPPE0030073.355		250 - 2500 µl	115 mm	Röd	1 box incl. 5 trays, 48 each	240

epT.I.P.S. box 2.

VWR art.nr	Antal	Tip style	Volym	Längd	Färgkod	Förp.	Förp.
EPPE0030076.125		Eppendorf Quality™	0,1 - 10 µl	34 mm	Mörkgrå	1 reusable box × 96 tips	96
EPPE0030076.133		Eppendorf Quality™	0,1 - 20 µl	40 mm	Grå	1 reusable box × 96 tips	96
EPPE0030076.214		Eppendorf Quality™	0,1 - 5000 µl	120 mm	Violett	1 reusable box × 24 tips	24
EPPE0030076.206		Eppendorf Quality™	0,25 - 2500 µl	115 mm	Röd	1 reusable box × 48 tips	48
EPPE0030076.141		Eppendorf Quality™	0,5 - 20 µl	46 mm	Ljusgrå	1 reusable box × 96 tips	96
EPPE0030076.150		Eppendorf Quality™	2 - 200 µl	53 mm	Gul	1 reusable box × 96 tips	96
EPPE0030076.168		Eppendorf Quality™	20 - 300 µl	55 mm	Orange	1 reusable box × 96 tips	96
EPPE0030076.176		Eppendorf Quality™	50 - 1000 µl	71 mm	Blå	1 reusable box × 96 tips	96
EPPE0030076.184		Eppendorf Quality™	50 - 1250 µl	76 mm	Grön	1 reusable box × 96 tips	96
EPPE0030076.192		Eppendorf Quality™	50 - 1250 µl	103 mm	Mörkgrön	1 reusable box × 96 tips	96

epT.I.P.S. set 2.

VWR art.nr	Antal	Tip style	Volym	Längd	Färgkod	Förp.	Förp.
EPPE0030076.290		Eppendorf Quality™	0,1 - 10 µl	34 mm	Mörkgrå	480 tips (5 trays × 96 tips), 1 reusable box	480
EPPE0030076.273		Eppendorf Quality™	0,1 - 20 µl	42 mm	Ljusrosa	1920 tips (5 trays × 384 tips), 1 reusable box	1.920
EPPE0030076.303		Eppendorf Quality™	0,1 - 20 µl	40 mm	Grå	480 tips (5 trays × 96 tips), 1 reusable box	480

EPPE0030076.362		Eppendorf Quality™	0,25 - 2500 µl	115 mm	Röd	240 tips (5 trays × 48 tips), 1 reusable box
EPPE0030076.311		Eppendorf Quality™	0,5 - 20 µl	46 mm	Ljusgrå	480 tips (5 trays × 96 tips), 1 reusable box
EPPE0030076.281		Eppendorf Quality™	0,5 - 100 µl	53 mm	Ljusgul	1920 tips (5 trays × 384 tips), 1 reusable box
EPPE0030076.320		Eppendorf Quality™	2 - 200 µl	53 mm	Gul	480 tips (5 trays × 96 tips), 1 reusable box
EPPE0030076.338		Eppendorf Quality™	20 - 300 µl	55 mm	Orange	480 tips (5 trays × 96 tips), 1 reusable box
EPPE0030076.346		Eppendorf Quality™	50 - 1000 µl	71 mm	Blå	480 tips (5 trays × 96 tips), 1 reusable box
EPPE0030076.354		Eppendorf Quality™	50 - 1250 µl	76 mm	Grön	480 tips (5 trays × 96 tips), 1 reusable box

Denna produkt är inte längre tillgänglig. Alternativ kan vara tillgängliga genom att söka med VWR art.nr ovan. Om du behöver ytterligare hjälp, ring VWR Kun