

CPR BARRIER DEVICES

The Laerdal® Pocket Mask

The Laerdal Pocket Mask's design, filter and one-way valve make it a quality choice that provides protection during CPR. By creating a slight distance between the patient and care provider the Pocket Mask also makes seeing patient lip color and chest movement easier.

A high-quality Technostat filter under a one-way valve filters out contaminants in the air being transferred between the care provider and the patient, protecting them both.



- 82000633 Pocket Mask in Re-sealable Polybag
- 82001133 Pocket Mask w/Gloves & Wipe in Hard Case
- 82001933 Pocket Mask only in Hard Case
- 82004033 Pocket Mask w/Gloves & Wipe in Blue Soft Case
- 82004133 Pocket Mask w/Gloves & Wipe in Black Soft Case
- 82004233 Pocket Mask w/Gloves & Wipe in Camo Soft Case
- 83001133 Pocket Mask with O2 Inlet, Headstrap, Gloves in Hard Case
- 83001633 Pocket Mask with O2 Inlet, Headstrap in Resealable Poly Bag
- 83001933 Pocket Mask with O2 Inlet, Headstrap in Hard Case
- 83004033 Pocket Mask with O2 Inlet, Headstrap, Gloves in Blue Soft Case
- 83004133 Pocket Mask with O2 Inlet, Headstrap, Gloves in Black Soft Case
- 83004233 Pocket Mask with O2 Inlet, Headstrap, Gloves in Camo Soft Case

The Laerdal Pocket Mask fits The BAG II and Laerdal Silicone Resuscitators.

[Click here for more information on the Laerdal Pocket Mask including directions for use.](#)

Consumables

- 820410 One-Way Valve (Qty 1)
- 820411 One-Way Valve for Pocket Mask (Qty 100)
- 820610 Filter Pocket Mask (Qty 1)
- 820611 Filter Pocket Mask (Qty 10)

