

# Pocket BVM

## Bag Valve Mask

---



---

Item No:	BR-00211
Manufacturer:	Micro BVM International Company Ltd.
Certificates:	CE
Classification:	N/A
Size:	21.7 x 12.1 cm
Color:	White
Country of origin:	CN
NSN no:	6515-01-590-8909
Tariff:	90192020
Weight:	450 g
Dimensions:	6.3 x 13.4 cm
Shelf life:	N/A

The Pocket BVM is an innovative resuscitator that provides breathing support to patients experiencing respiratory distress. Unlike any other BVM, the Pocket BVM has a patented collapsible design that enables it to fit into a compact and durable case. This unique design makes it the smallest fully functional BVM on the market.

Since the Pocket BVM is so small (13.4 cm diameter x 6.3 cm height) it easily fits into any emergency kit. This is an especially important factor for paramedics that have space issues such as those on bikes, motorcycles, all-terrain vehicles (ATVs), helicopters etc. The Pocket BVM also saves valuable space, and on average uses a quarter of the space compared to standard BVMs.