

KV-40 1/2 CIRCLE, 48 mm



Order Information

Code	Trademark	Length	Color	Box Qty	U.S.P. (E.P.)	Strands Per Pkg
CL840	POLYSORB™	30in (75cm)	VIOLET	36	1 (4)	1
CL841	POLYSORB™	30in (75cm)	VIOLET	36	2 (5)	1
CL941	POLYSORB™	36in (90cm)	VIOLET	36	2 (5)	1
8886239369	STEEL	18in (45cm)	NONE	12	6 (8)	4
8886239389	STEEL	18in (45cm)	NONE	12	5 (7)	4
8886239269	STEEL	18in (45cm)	NONE	12	6 (8)	2
8886239249	STEEL	18in (45cm)	NONE	12	7 (9)	2
8886241489	STEEL	18in (45cm)	NONE	12	5 (7)	4
8886281889	TI•CRON™	30in (75cm)	BLUE	12	5 (7)	4
8886246625	STEEL	18in (45cm)	NONE	12	3 (6)	5
8886629571	MAXON™	36in (90cm)	GREEN	36	1 (4)	1
CL739	POLYSORB™	36in (90cm)	VIOLET	36	0 (3.5)	1
CL940	POLYSORB™	36in (90cm)	VIOLET	36	1 (4)	1
CL740L	POLYSORB™	60in (150cm)	VIOLET	24	1 (4)	1
CL741L	POLYSORB™	60in (150cm)	VIOLET	24	2 (5)	1
XX5097L	TI•CRON™	60in (150cm)	BLUE	36	1 (4)	1
CM740L	BIOSYN™	60in (150cm)	VIOLET	12	1 (4)	1
GMM807L	MAXON™	60in (150cm)	GREEN	12	1 (4)	1
GMM806L	MAXON™	60in (150cm)	GREEN	12	0 (3.5)	1