


C-14 3/8 CIRCLE, 24 mm



Order Information

Code	Trademark	Double Arm	Length	Color	Box Qty	U.S.P. (E.P.)	Strands Per Pkg
8886193433	SURGILON™	No	18in (45cm)	BLACK	12	4-0 (1.5)	1
8886193443	SURGILON™	No	18in (45cm)	BLACK	12	3-0 (2)	1
8886193453	SURGILON™	No	18in (45cm)	BLACK	12	2-0 (3)	1
CL823	POLYSORB™	No	30in (75cm)	VIOLET	36	2-0 (3)	1
8886305041	TI•CRON™	No	30in (75cm)	BLUE	36	3-0 (2)	1
SG637	CHROMIC GUT	No	30in (75cm)	NONE	36	2-0 (3.5)	1
SG637G	CHROMIC GUT	No	30in (75cm)	NONE	12	2-0 (3.5)	1
SG810	PLAIN GUT	No	30in (75cm)	NONE	36	3-0 (3)	1
SG823	PLAIN GUT	No	30in (75cm)	NONE	36	2-0 (3.5)	1
8886174431	DERMALON™	No	18in (45cm)	BLUE	36	4-0 (1.5)	1
8886179941	DERMALON™	No	30in (75cm)	BLUE	36	3-0 (2)	1
8886174441	DERMALON™	No	18in (45cm)	BLUE	36	3-0 (2)	1
SN663G	MONOSOF™	No	18in (45cm)	BLACK	12	3-0 (2)	1
SN663	MONOSOF™	No	18in (45cm)	BLACK	36	3-0 (2)	1
SN669G	MONOSOF™	No	30in (75cm)	BLACK	12	3-0 (2)	1
SN669	MONOSOF™	No	30in (75cm)	BLACK	36	3-0 (2)	1
SN658	MONOSOF™	No	18in (45cm)	BLACK	36	4-0 (1.5)	1
SN649	MONOSOF™	No	36in (90cm)	BLACK	36	4-0 (1.5)	1
SN648	MONOSOF™	No	36in (90cm)	BLACK	36	3-0 (2)	1
8886442253	NOVAFIL™	No	18in (45cm)	BLUE	12	2-0 (3)	1
8886442243	NOVAFIL™	No	18in (45cm)	BLUE	12	3-0 (2)	1
8886442233	NOVAFIL™	No	18in (45cm)	BLUE	12	4-0 (1.5)	1
SP629	SURGIPRO™	No	18in (45cm)	BLUE	36	4-0 (1.5)	1
SP684	SURGIPRO™ II	No	18in (45cm)	BLUE	36	3-0 (2)	1
SP684G	SURGIPRO™	No	18in (45cm)	BLUE	12	3-0 (2)	1
SS684	SOFSILK™	No	18in (45cm)	BLACK	36	3-0 (2)	1
SS629	SOFSILK™	No	18in (45cm)	BLACK	36	4-0 (1.5)	1
SS684G	SOFSILK™	No	18in (45cm)	BLACK	12	3-0 (2)	1
SL692G	POLYSORB™	No	30in (75cm)	UNDYED	12	4-0 (1.5)	1
SL694	POLYSORB™	No	30in (75cm)	VIOLET	36	3-0 (2)	1
SL823	POLYSORB™	No	30in (75cm)	UNDYED	36	2-0 (3)	1
SL693	POLYSORB™	No	30in (75cm)	UNDYED	36	3-0 (2)	1
SL692	POLYSORB™	No	30in (75cm)	UNDYED	36	4-0 (1.5)	1
SM823G	BIOSYN™	No	30in (75cm)	UNDYED	12	2-0 (3)	1
SM823	BIOSYN™	No	30in (75cm)	UNDYED	36	2-0 (3)	1
SM693G	BIOSYN™	No	30in (75cm)	UNDYED	12	3-0 (2)	1

SM693	BIOSYN™	No	30in (75cm)	UNDYED	36	3-0 (2)	1
8886664151	MAXON™	No	30in (75cm)	GREEN	36	2-0 (3)	1
8886664141	MAXON™	No	30in (75cm)	GREEN	36	3-0 (2)	1
8886661851	MAXON™	No	30in (75cm)	CLEAR	36	2-0 (3)	1
8886661841	MAXON™	No	30in (75cm)	CLEAR	36	3-0 (2)	1
8886661831	MAXON™	No	30in (75cm)	CLEAR	36	4-0 (1.5)	1
8886306561	TI•CRON™	No	36in (90cm)	BLUE	36	0 (3.5)	1
8886306941	TI•CRON™	No	18in (45cm)	BLUE	36	3-0 (2)	1
8886193541	SURGILON™	No	30in (75cm)	BLACK	36	3-0 (2)	1
8886193551	SURGILON™	No	30in (75cm)	BLACK	36	2-0 (3)	1
SN758	MONOSOF™	No	18in (45cm)	BLACK	36	2-0 (3)	1
SC693G	CAPROSYN™	No	30in (75cm)	UNDYED	12	3-0 (2)	1
8886306551	TI•CRON™	No	36in (90cm)	BLUE	36	2-0 (3)	1
SS695	SOFSILK™	No	18in (45cm)	BLACK	36	2-0 (3)	1
8886442153	NOVAFIL™	No	30in (75cm)	BLUE	12	2-0 (3)	1
8886305051	TI•CRON™	No	30in (75cm)	BLUE	36	2-0 (3)	1
SS784	SOFSILK™	No	30in (75cm)	BLACK	36	3-0 (2)	1
SS783	SOFSILK™	No	30in (75cm)	BLACK	36	4-0 (1.5)	1
SS786	SOFSILK™	No	30in (75cm)	BLACK	36	0 (3.5)	1
SS785	SOFSILK™	No	30in (75cm)	BLACK	36	2-0 (3)	1
8886193561	SURGILON™	No	30in (75cm)	BLACK	36	0 (3.5)	1
SN747	MONOSOF™	No	36in (90cm)	BLACK	36	2-0 (3)	1
SS787	SOFSILK™	No	30in (75cm)	BLACK	36	1 (4)	1
8886442263	NOVAFIL™	No	18in (45cm)	BLUE	12	0 (3.5)	1
8886400241	NOVAFIL™	No	36in (90cm)	BLUE	36	3-0 (2)	1
X1079	BIOSYN™	No	30in (75cm)	VIOLET	36	2-0 (3)	1
8886907551	DEXON™	No	30in (75cm)	BEIGE	36	2-0 (3)	1
SM923	BIOSYN™	No	30in (75cm)	VIOLET	36	2-0 (3)	1
SL695	POLYSORB™	No	30in (75cm)	VIOLET	36	4-0 (1.5)	1
CL900	POLYSORB™	GS-21	60in (150cm)	VIOLET	36	2-0 (3)	1
8886305061	TI•CRON™	No	30in (75cm)	BLUE	36	0 (3.5)	1
8886179951	DERMALON™	No	30in (75cm)	BLUE	36	2-0 (3)	1
8886179931	DERMALON™	No	30in (75cm)	BLUE	36	4-0 (1.5)	1
SM1023G	BIOSYN™	No	36in (90cm)	VIOLET	12	2-0 (3)	1