





**Polysorb™ Lactomer™ 9-1, Coated, Braided Synthetic Absorbable Sutures**

Reverse Cutting								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
<b>1/2 CIRCLE</b> 40 mm H 	2-0	3	5x45	Violet	12	D-Tach™		CL82MG
	0	3.5	5x45	Undyed	12	D-Tach™		CL76MG
	0	3.5	5x45	Violet	12	D-Tach™		CL80MG
	1	4	5x45	Undyed	12	D-Tach™		CL77MG
	1	4	5x45	Violet	12	D-Tach™		CL81MG
<b>1/2 CIRCLE</b> 48 mm 	2-0	3	75	Violet	36			CL571
	0	3.5	75	Violet	36			CL568
	0	3.5	90	Violet	36			CL927
	1	4	75	Undyed	36			CL555
	1	4	75	Violet	36			CL569
	1	4	90	Violet	36			CL928
	2	5	75	Undyed	36			CL960
	2	5	75	Violet	36			CL570
	2	5	90	Undyed	36			CL939
	2	5	90	Violet	36			CL929
<b>1/2 CERCHIO</b> 55 mm Super Rinforzato 	2	5	75	Violet	36			CL979
<b>STRAIGHT</b> 60 mm 	4-0	1.5	75	Undyed	36			SL644
	4-0	1.5	75	Violet	36			SL654
	3-0	2	75	Undyed	36			SL643
	3-0	2	75	Violet	36			SL653
	2-0	3	75	Undyed	36			SL642
	2-0	3	75	Violet	36			SL652
	2-0	3	75	Violet	12			SL652G