






Taper Point								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code	
1/2 CIRCLE 26 mm 	0	3.5	75	Violet	36			GL124
	1	4	75	Violet	36			GL130
1/2 CIRCLE 26 mm   STRAIGHT 44 mm 	3-0	2	90	Undyed	36			GL423
1/2 CIRCLE 26 mm 	3-0	2	90	Violet	36			GL332
	2-0	3	90	Violet	36			GL333
1/2 CIRCLE 27 mm 	3-0	2	75	Undyed	36			CL829
	3-0	2	75	Violet	36			CL882
	2-0	3	5x45	Violet	12	D-Tach™		CL16MG
	2-0	3	75	Undyed	36	Surgalloy™		CL5830
	2-0	3	75	Undyed	36			CL830
	2-0	3	75	Violet	36	Surgalloy™		CL5883
	2-0	3	75	Violet	36			CL883
	0	3.5	5x45	Violet	12	D-Tach™		CL17MG
	0	3.5	75	Undyed	36	Surgalloy™		CL5831
	0	3.5	75	Undyed	36			CL831
	0	3.5	75	Violet	36	Surgalloy™		CL5884
	0	3.5	75	Violet	36			CL884
	0	3.5	5x75	Violet	12	D-Tach™		CL886M