




Polysorb™ Lactomer™ 9-1, Coated, Braided Synthetic Absorbable Sutures

Taper Point							
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Plating-out	Reorder code
1/2 CIRCLE 22 mm 	3-0	2	90	Violet	36		GL282
1/2 CIRCLE 22 mm 	3-0	2	75	Undyed	36		CL896
	3-0	2	75	Violet	36		CL892
	2-0	3	75	Violet	36		CL893
	0	3.5	75	Violet	36		CL894
1/2 CIRCLE 26 mm 	5-0	1	75	Violet	36		GL120
	4-0	1.5	5x45	Violet	12	D-Tach™	GL64MG
	4-0	1.5	75	Undyed	36	Surgalloy™	GL5321
	4-0	1.5	75	Undyed	36		GL321
	4-0	1.5	75	Violet	36	Surgalloy™	GL5121
	4-0	1.5	75	Violet	36		GL121
	4-0	1.5	5x75	Violet	12	D-Tach™	GL68M
	3-0	2	5x45	Undyed	12	D-Tach™	GL69MG
	3-0	2	5x45	Violet	12	D-Tach™	GL63MG
	3-0	2	75	Undyed	36	Surgalloy™	GL5322
	3-0	2	75	Undyed	36		GL322
	3-0	2	75	Violet	36	Surgalloy™	GL5122
	3-0	2	75	Violet	36		GL122
	3-0	2	5x75	Violet	12	D-Tach™	GL67M
	3-0	2	90	Violet	36		GL226
	2-0	3	5x45	Undyed	12	D-Tach™	GL61MG
	2-0	3	5x45	Violet	12	D-Tach™	GL62MG
	2-0	3	75	Undyed	36	Surgalloy™	GL5323
	2-0	3	75	Undyed	36		GL323
	2-0	3	75	Violet	36	Surgalloy™	GL5123
	2-0	3	75	Violet	36		GL123
	2-0	3	5x75	Violet	12	D-Tach™	GL66M
	2-0	3	90	Violet	36		GL125
	0	3.5	75	Undyed	36		GL324