

Shiley™ Oral/Nasal Tracheal Tube, Cuffless



Features and benefits

- Smooth beveled and carefully molded hooded tip
- Graduation marks on both sides of the tracheal tube and special markings at the tip
- Tip-to-tip x-ray line
- Incorporated Murphy eye
- Suitable for both oral and nasal intubation

Ordering information

	REF	I.D. (mm)	O.D. (mm)	LENGTH (mm)
PEDIATRIC	86232	2.0	2.9	171
	86233	2.5	3.6	171
	86234	3.0	4.3	190
	86235	3.5	4.9	213
	86236	4.0	5.6	232
	86237	4.5	6.2	253
	86238	5.0	6.9	274
	86239	5.5	7.5	304
	86240	6.0	8.2	317
	86241	6.5	8.8	327
86242	7.0	9.6	339	

Carton quantity: 10

All tubes are delivered sterile packed, for single use only.

Shiley™ Contour Oral/Nasal Tracheal Tube, Cuffless

[LEARN MORE](#)



Features and benefits

- Soft material — May be suitable for nasal intubation
- Semitransparent tubes for direct visualization of exhaled air condensation
- Specially molded hooded tip
- Graduation marks on both sides of the tracheal tube and special markings at the tip
- Tip-to-tip x-ray line
- Incorporated Murphy eye

Ordering information

	REF	I.D. (mm)	O.D. (mm)	LENGTH (mm)
PEDIATRIC	411-25	2.5	3.6	171
	411-30	3.0	4.3	190
	411-35	3.5	5.1	213
	411-40	4.0	6.0	232
	411-45	4.5	6.5	253
	411-50	5.0	6.9	274
	411-55	5.5	7.5	304
	411-60	6.0	8.2	317
	411-65	6.5	8.8	327
	411-70	7.0	9.6	339

Carton quantity: 10

All tubes are delivered sterile packed, for single use only.