


<b>Product Name</b>	Straight Connector 22M-22M	
<b>Reference numbers</b>	1960000	
<b>Product Type</b>	Connector	
<b>Manufacturer</b>	Intersurgical	
<b>Materials</b>	PP	
<b>Packaging</b>	Individually packaged	
<b>Unit Packaging</b>	40	
<b>Box quantity</b>	10240	
<b>Pallet quantity</b>		
<b>Instructions for Use</b>	Instructions for use with product	N/A
	Instructions for use with each box	N/A
<b>Certification</b>	CE0120	Class 2A
<b>Applicable Standards</b>	BS EN ISO 13485:2012 Medical devices — Quality management systems — Requirements for regulatory purposes	
	ISO 9001:2008 Quality Management Systems - Requirements	
	ISO 15223-1:2012 Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied	
	BS EN ISO 5356-1:2015 - Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets	
<b>Sterile</b>	No	
<b>Shelf life</b>	5 years	
<b>Disposal Methods</b>	As hospital guidelines	
<b>Re-use</b>	ⓧ	Single use only

**Product characteristics**

Technical Information	
Weight (g)	6.2
Internal volume (ml)	13.9



A 3D perspective view of a green, cylindrical straight connector. The connector has a central lumen and a flange-like structure in the middle. Both ends are labeled '22M' with double-headed arrows indicating the connection size.

