

# Petri Dishes, Inoculation Loops



## Petri dishes for bacteriology

Made from clear polystyrene, heat resistant up to 80° C.

Order Number	Ø / Height mm	Ventilation Cams	Packaging Bag/Case
82.1184	150 / 20 gamma-sterile	with	10 / 100
82.1472	92 / 16	without	20 / 480
82.1473	92 / 16	with	20 / 480
82.1472.001	92 / 16 gamma-sterile	without	20 / 480
82.1473.001	92 / 16 gamma-sterile	with	20 / 480
82.1194	60 / 15 gamma-sterile	with	20 / 500
82.1135	35 / 10 gamma-sterile	with	20 / 500

## NEW Petri dishes, coloured, with ventilation cams

Order Number	Ø / Height mm	Colour	Packaging Bag/Case
82.1473.020	92 / 16	red	20 / 480
82.1473.040	92 / 16	yellow	20 / 480
82.1473.060	92 / 16	blue	20 / 480
82.1473.080	92 / 16	green	20 / 480

## Petri dish, round, with 2 compartments

Order Number	Ø / Height mm	Ventilation Cams	Packaging Bag/Case
82.1195	92 / 16	with	20 / 480

## Petri dish, square


Order Number	Length x Width x Height mm	Ventilation Cams	Packaging Bag/Case
82.9923.422	100 x 100 x 20 gamma-sterile	without	4 / 160

DishRack for Petri dishes 92/16 on p. 117.

## Gamma-sterile inoculation loops, inoculation needle and inoculation spreader from polystyrene

- Time-saving for laboratories with high work loads
- Easy to use
- Safe – no pathogenic aerosol and cross contamination

Order Number	Description	Colour	Packaging
86.1562.010	10 µl loop	blue	10/Bag, sterile 1,000/Case
86.1562.050	10 µl loop	blue	50/Bag, sterile 2,000/Case
86.1567.010	1 µl loop	white	10/Bag, sterile 1,000/Case
86.1567.050	1 µl loop	white	50/Bag, sterile 2,000/Case
86.1568.010	needle	orange	10/Bag, sterile 1,000/Case
86.1568.050	needle	orange	50/Bag, sterile 2,000/Case
86.1569.001	inoculation spreader	blue	1/Bag, sterile 500/Case
86.1569.005	inoculation spreader	blue	5/Bag, sterile 500/Case



NEW

NEW