

PRODUKTINFORMATION KANYL TILL SUBCUTAN VENPORT, GRIPPER®



KANYL TILL SUBCUTAN VENPORT, GRIPPER®

- Lätt att fixera
- Ett avtagbart grepp säkerställer en kontrollerad nålplacering
- En icke absorberande skumdyna för ökad patientkomfort
- Integrerad anslutnings slang ca 20 cm med luerkoppling och färgkodad avstängningsklämma
- Kanylslipning som minimerar risk för urstansning av portmembranet
- Ger en bra patientkomfort vid avlägsnande.
- Rekommenderad användningstid enl amerikansk standard från INS är 7 dagar
- Icke borrarande nål som saknar silikonbeläggning

- » CE märkt
- » Steriliseringsmetod EO
- » Engångsmodell
- » Innehåller ej Latex

Benämning	Vårt Art nr	Enhet	Avdfp	Trpf
Gripper®nål 19Gx19mm	21-2735-24	st	12	12
Gripper®nål 19Gx25mm	21-2720-24	st	12	12
Gripper®nål 19Gx32mm	21-2721-24	st	12	12
Gripper®nål 20Gx16mm	21-2736-24	st	12	12
Gripper®nål 20Gx19mm	21-2734-24	st	12	12
Gripper®nål 20Gx25mm	21-2717-24	st	12	12
Gripper®nål 20Gx32mm	21-2718-24	st	12	12
Gripper®nål 22Gx16mm	21-2737-24	st	12	12
Gripper®nål 22Gx19mm	21-2733-24	st	12	12
Gripper®nål 22Gx25mm	21-2714-24	st	12	12
Gripper®nål 22Gx32mm	21-2715-24	st	12	12

Smiths Medical Sverige AB.
Orrvägen 26-28, 192 55 Sollentuna.
Tel. växel 08- 594 772 50
Fax. 08-594 772 69
www.smiths-medical.com

Vid frågor kontakta er
lokala säljrepresentant

Smiths Medical is part of the global technology business Smiths Group plc. Product(s) described may not be licensed or available for sale in all countries. GRIPPER PLUS, Deltec and the Smiths Medical design marks are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. ©2012 Smiths Medical. All rights reserved. XX19XXXXSV-XXXX