

PRODUKTINFORMATION - TUOHYNÅL



TUOHYNÅL MED TILLHÖRANDE VINGE

- Optimerad Tuohy slipning ger maximal känsla vid penetration
- Distansmarkeringar på 10mm
- För att minska risken för "vävnadsproppar" i Tuohynålen är mandrinen individuellt anpassad till varje enskild Tuohynål
- Färgkodad epidural Tuohynål för säker och enkel identifiering; Ljusblå färgkod är 16G, mörkblå färgkod är 18G
- Har anslutning av luerfattning
- Löstagbar vinge för ökad flexibilitet.

» Latexfri

» CE märkt

» Steriliseringsmetod E0

» Engångsartikel

Benämning	Vårt Art nr	Gauge	Längd	Förpackning
Tuohynål med separat tillhörande vinge, 16G, ljusblå	100/395/560	16G	110mm	10
Tuohynål med separat tillhörande vinge, 16G, ljusblå	100/395/580	18G	110mm	10

PORTEX®

Smiths Medical Sverige AB
Box 1143, 164 22 Kista
Tel. 08- 594 772 50
Fax. 08-594 772 69
www.smiths-medical.com

Vid frågor kontakta er
lokala säljrepresentant

Smiths Medical is part of the global technology business Smiths Group plc. Product(s) described may not be licensed or available for sale in all countries. LEVEL 1 and the Smiths Medical design marks are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. ©2012 Smiths Medical. All rights reserved.
20120917SV-0001

smiths medical