




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/uk



Technical sheet

Surgical gown

Reference(s)

9505CEA, 9506CEA, 9507CEA, 9516CEA, 9575CEA

Description

A sterile wrap 60x60 cm contains:

- 1 Optima® surgical gown, Velcro® adjustable neck
- 2 hand towels, 39x30 cm unfolded

Reference	Length (cm)	Size	Unit/case
9505CEA	115	S/M	28
9506CEA	120	L	28
9507CEA	130	XL	24
9516CEA	150	XLL	24
9575CEA	133	XXL	24

Main materials

Optima® 72g/m² (polyester + woodpulp)

Packaging

Sterilization pouch

2 polyliners contain all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile



Latex free



0123



STERILE R

© 2011 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-95000CE-EN01

01.05.SUB.04