

Discofix® C – Settle for the best

I.V. Accessories & Filters



Product Type		Version		Sales unit/Pcs.	Code No. (REF)
C-3	blue	Single sterile		100	16494C
C-3	white	Single sterile		100	16495C
C-3	red	Single sterile		100	16496C
C-3	yellow	Single sterile		100	16497C
C-3	blue	Single sterile with extension	10 cm	50	16500C
C-3	blue	Single sterile with extension	25 cm	50	16520C
C-3	blue	Single sterile with extension	50 cm	50	16540C
C-3	blue	Single sterile with extension	75 cm	50	16551C
C-3	blue	Single sterile with extension	100 cm	50	16560C
C-3	white	Single sterile with extension	10 cm	50	16501C
PinPad		Easy clip-on fixation system, fits to Discofix C		25	4515510
C-3 Manifold	b/y/r *	Single sterile		50	16600C
C-3 Manifold	blue	Single sterile		50	16605C
C-3 Manifold	blue	Single sterile with extension and proximal blue stopcock	150 cm	50	16610C
C-3 Manifold	blue	Single sterile with extension	150 cm	50	16615C
		Industrial customers		Sales unit/Pcs.	Code No. (REF)
C-3	blue	Bulk non-sterile		on request	16016C
C-3	red	Bulk non-sterile		on request	16017C
C-3	white	Bulk non-sterile		on request	16018C
C-3 Manifold	b/y/r *	Bulk non-sterile		on request	16620C
C-3 Manifold	blue	Bulk non-sterile		on request	16625C

* blue/yellow/red



B. Braun Melsungen AG

P.O. Box 1120
D-34209 Melsungen
Tel +49(0) 56 61 71-0

www.bbraun.com