

Product Data Sheet

Spruta luerlock 3komp Terumo

REF: SS*10LE1



Tryckt: 29.4.2016

REF	Storlek	
SS*10LE1	10 mL	
CE 0197		A photograph showing five white plastic syringes of different sizes, arranged in a descending diagonal line from top-left to bottom-right. Each syringe has a blue plunger and a clear barrel with markings. The largest syringe is at the top, and the smallest is at the bottom right.
Material:		Polypropylen, termoplastisk elastomer, silikonolja. Latex och PVC-fri.
Egenskaper:		Steril genomskinlig 3-komponent spruta med 0,2 ml
Steril produkt:		Ja
Hållbarhetstid:		5 år
MBA/st:		100
Avdelningsförpackning/st:		100
Transportlåda/st:		800