

## CORR-A-FLEX II Roll Tubing

- For heated and non-heated therapy applications
- Slips on and locks securely to 22-mm adaptors, thanks to collar on tubing ends
- Extremely lightweight and flexible
- With cuttable sections every 15 cm
- 30.4-m roll; polyurethane (PE)



**CORR-A-FLEX II Roll Tubing**

Reference	Description	Pack factor
1680	Roll tubing	1 roll

## Aerosol Drainage System

- Self-locking plastic chain secures bag to a stable object
- Disposable drainage bag with a tee adaptor which accepts standard 22-mm ID corrugated tubing
- Convenient drainage port
- Material: Tee polyethylene, bag PVC




**Aerosol Drainage System**

Reference	Description	Pack factor
41740	750 ml	50

## CORR-A-FLEX tube

- Collared tubing end slips on and locks to 22 mm adaptors for secure fit
- Exceptionally flexible, lightweight tubing, PE
- Length of 180 cm
- With cuttable sections every 15 cm
- Individually packed



**CORR-A-FLEX Tube**

Reference	Description	Pack factor
41518	Roll tubing	50

## ADDIPAK Unit Dose

- 3-ml and 15-ml vials of sterile 0.9% saline solution for inhalation therapy
- Inverted millilitre graduations for accurate dispensing



**ADDIPAK Unit Dose**

Reference	Description	Pack factor
200-39	0.9% - full normal saline, 3 ml (pink vial)	100
200-69	0.9% - full normal saline, 15 ml (pink vial)	48