

SKINTACT Diagnostic Resting ECG Electrodes

"SKINTACT Easitabs are designed for Resting 12 lead diagnostic ECG's. Easitabs are disposable eliminating the spread of cross infection associated with older style suction cups. As the name suggests Easitabs are easy to use, provide reliable high quality recordings, and are suitable for all types of patients."



Manufactured using Leonhard Lang's unique C-LINE technology which ensures that every electrode will be fresh providing you with exceptional consistency in quality of recordings at very affordable prices.

SKINTACT Diagnostic Resting ECG Electrodes

The SKINTACT Easitab System –

Your best choice for an excellent and safe resting ECG:

Excellent Trace Quality. The conductive adhesive Aqua-Tac hydrogel of SKINTACT Easitabs provides a reliable signal tracing of consistent high quality. Easitab FS-WB00 completes the range for demanding conditions.

Secure Adhesion and Patient Comfort. SKINTACT Easitabs feel soft but stick well to the patient's skin. Their contact design prevents lift-off. After the ECG they can be removed comfortably – without loss of hair.

Superior Safety. Proper use of SKINTACT Easitabs eliminates cross infection whereas suction electrode bred bacteria and fungi.


Time Efficiency. Unlike suction electrodes, SKINTACT Easitabs on a handy card are applied quickly on all types of patients. Aqua-Tac gel requires no clean up after the ECG. Easiprep and Easiclip tab adaptors further facilitate your success.

As an alternative we offer you a choice of SKINTACT F-RG with standard stud connector and SKINTACT FS-RB4 with adaptor for a 4 mm banana connector *).



Order Code	RT-14	RT-34	RT-74	FS-WB00	F-RG	FS-RB4 *)
Gel, Connector	Aqua-Tac, Easitab	Aqua-Tac, Easitab	Aqua-Tac, Easitab	Aqua-Wet, Easitab	Aqua-Wet, Snap	Aqua-Wet, 4 mm Adaptor
Backing, Adhesive	Paper/PET, Gel	Paper/PET, Gel	Paper/PET, Gel	Foam, Enhanced	Foam, Standard	Foam, Enhanced
Dimensions	14 x 34 mm pediatric, card of 10	22 x 34 mm, adult, card of 10	26 x 34 mm, adult, card of 10	28 x 37 mm, adult, card of 10	32 x 41 mm, universal, card of 6	28 x 44 mm, adult, card of 3
Packaging Standard	100/pouch 500/box 5.000/case	100/pouch 500/box 5.000/case	100/pouch 500/box 5.000/case	50/pouch 500/box 3.000/case	30/pouch 1.200/box 4.800/case	30/pouch 750/box 3.000/case

*) Also available with 3 mm adaptor (order code FS-RB3).



Product	Easiclip	Easiprep	Easiprep	Easiclean
Order Code	58100	58727	58726	59635
Packaging Standard	10/pouch	500/pouch	50/pouch	100/box

All SKINTACT diagnostic resting ECG electrodes meet the ANSI/AAMI standard EC12:2000 for disposable pregelled ECG electrodes.