

Disposable Containers

with snap-on lid



250 ml Container, polystyrene, white

Order Number	Length mm	Aperture-Ø mm	Packaging Pcs./Case
75.560	100	73	1,550

Lid for 250 ml container

Order Number: 76.561

Packaging: 1,550/Case



125 ml Container, neutral polypropylene

Order Number	Length mm	Aperture-Ø mm	Packaging Pcs./Case
75.565	85	62	500

Lid for 125 ml container

Order Number: 76.566

Packaging: 2,000/Case



125 ml Container, polypropylene, highly transparent

Order Number	Length mm	Aperture-Ø mm	Packaging Pcs./Case
75.1355	67	66	1,000
75.1355.100	67	66	100/Bag / 1,000/Case

Snap-on lid, red polyethylene, for 125 ml container

Order Number	Description	Pcs./Bag	Packaging Pcs./Case
76.1477	with spout	1,000	1,000
76.1477.100	with spout	100	1,000
76.1414	without spout	1,000	1,000
76.1414.100	without spout	100	1,000




100 ml Container, polystyrene, white

Order Number	Length mm	Aperture-Ø mm	Packaging Pcs./Case
75.570	49	70	1,000

Snap-on lid for 100 ml container

Order Number: 76.571

Packaging: 1,000/Case



30 ml Medicine cup from polypropylene, with graduation

Order Number	Length/Aperture-Ø mm	Colour	Packaging Pcs./Case
75.1337	40/37	neutral	3,600 each case (60 per Bag)
75.1337.002	40/37	red	
75.1337.004	40/37	yellow	
75.1337.005	40/37	green	
75.1337.006	40/37	blue	

NEW
NEW
NEW
NEW

Snap-on lid for 30 ml medicine cup, white polyethylene

Order Number: 76.1340

3,600/Case

Order Number: 76.1340.060

60/Bag, 3,600/Case

PE = polyethylene · PP = polypropylene · PS = polystyrene