

3M

Littmann[™] Stethoscopes

Parts & Accessories Reference Guide




Eartips

3M™ Littmann™ eartips have been improved to offer increased comfort, durability and a tighter fit to your Littmann stethoscope.

Ribbed eartubes on all Littmann stethoscopes manufactured after May 1994 offer a “Snap-tight” fit that will not loosen during routine use. Littmann stethoscopes manufactured before May 1994 feature a threaded eartip design.

“Threaded” or “Snap-tight” eartips?
How to determine which eartips to order:


Diagram A



Order “Threaded” eartips

If your stethoscope has threaded eartubes similar to the adjacent diagram A
-or-
If the eartips being replaced have a plastic threaded insert

Diagram B



Order “Snap-tight” eartips

If your stethoscope has ribbed eartubes similar to the adjacent diagram B
-or-
If the eartips being replaced have no plastic threaded insert and an end opening of 1/8" diameter

3M™ Littmann™ Eartips

Description	Color	Size	Order Qty.	Cat. No.	11 Digit Number
-------------	-------	------	------------	----------	-----------------

Soft-sealing eartips – snap-tight / push into place (for Littmann stethoscopes manufactured after May 1994)

1	Soft-sealing eartips Snap-tight	Black	Large	10 pair/carton	37811	78-8078-3525-7
	Soft-sealing eartips Snap-tight	Gray	Large	10 pair/carton	37812	78-8078-3526-5
2	Soft-sealing eartips Snap-tight	Black	Small	10 pair/carton	37809	78-8078-3527-3
	Soft-sealing eartips Snap-tight	Gray	Small	10 pair/carton	37810	78-8078-3528-1

Soft-sealing eartips – threaded / screw into place (for Littmann stethoscopes manufactured before May 1994)

1	Soft-sealing eartips Threaded	Black	Large	10 pair/carton	36534	78-8078-2775-9
	Soft-sealing eartips Threaded	Gray	Large	10 pair/carton	36536	78-8078-2777-5
2	Soft-sealing eartips Threaded	Black	Small	10 pair/carton	36535	78-8078-2776-7
	Soft-sealing eartips Threaded	Gray	Small	10 pair/carton	36537	78-8078-2778-3

Firm eartips – snap-tight / push into place (for Littmann stethoscopes manufactured after May 1994)

3	Firm eartips Snap-tight	Gray	Large	10 pair/bag	37808	78-8078-3529-9
4	Firm eartips Snap-tight	Gray	Small	10 pair/bag	38742	78-8078-4103-2

Firm eartips – threaded / screw into place (for Littmann stethoscopes manufactured before May 1994)


3	Firm eartips Threaded	Gray	Large	10 pair/bag	36538	78-8078-2779-1
4	Firm eartips Threaded	Gray	Small	10 pair/bag	36539	78-8078-2780-9

Hard eartips – threaded / screw into place (for Littmann stethoscopes manufactured before May 1994)


5	Hard eartips Threaded	White	Large	5 pair/bag	36540	78-8078-2781-7
6	Hard eartips Threaded	White	Small	5 pair/bag	36541	78-8078-2782-5

Eartips (cont.)


Soft-sealing eartips







Firm eartips



Hard eartips



Nonchill Bell Sleeves




Description	Color	Order Qty.	Cat. No.	11 Digit Number
-------------	-------	------------	----------	-----------------

Nonchill bell sleeve for

7	Cardiology III / Cardiology II S.E. / Cardiology II	Black	10/bag	36570	78-8078-7338-1
	Cardiology III / Cardiology II S.E. / Cardiology II	Gray	10/bag	36571	78-8078-7337-3
8	Classic II S.E. / Classic II / Lightweight	Black	10/bag	36544	78-8078-2785-8
	Classic II S.E. / Classic II / Lightweight	Gray	10/bag	36545	78-8078-2786-6
9	Classic II Pediatric	Black	5/bag	36546	78-8078-2787-4
	Classic II Pediatric	Gray	5/bag	36568	78-8078-8013-9
10	Classic II Infant	Black	5/bag	36547	78-8078-2788-2
	Classic II Infant	Gray	5/bag	36569	78-8078-8014-7

Diaphragms and Rim Assemblies



Description	Color	Order Qty.	Cat. No.	11 Digit Number
-------------	-------	------------	----------	-----------------

Tunable Diaphragm & Rim Assembly for

1	Master Cardiology	Black Rim	5/bag	36554	78-8078-2795-7
	Master Cardiology	Gray Rim	5/bag	36555	78-8078-2796-5
2	Cardiology III (large side)	Black Rim	5/bag	36556	78-8078-2797-3
	Cardiology III (large side)	Gray Rim	5/bag	36557	78-8078-2798-1
3	Cardiology III (small side)	Black Rim	5/bag	36572	78-8078-7340-7
	Cardiology III (small side)	Gray Rim	5/bag	36573	78-8078-7339-9
2	Cardiology II S.E.	Black Rim	5/bag	36556	78-8078-2797-3
	Cardiology II S.E.	Gray Rim	5/bag	36557	78-8078-2798-1
2	Master Classic II	Black Rim	5/bag	36556	78-8078-2797-3
	Master Classic II	Gray Rim	5/bag	36557	78-8078-2798-1
2	Classic II S.E.	Black Rim	5/bag	36556	78-8078-2797-3
	Classic II S.E.	Gray Rim	5/bag	36557	78-8078-2798-1
2	Select	Black Rim	5/bag	36556	78-8078-2797-3
	Select	Gray Rim	5/bag	36557	78-8078-2798-1

Diaphragm (Printed with the Littmann Logo) for


4	Cardiology II		10/bag	36559	78-8078-2799-9
4	Classic II		10/bag	36559	78-8078-2799-9
4	Lightweight		10/bag	36559	78-8078-2799-9
5	Classic II Pediatric		10/bag	36560	78-8078-2800-5
6	Classic II Infant		10/bag	36561	78-8078-2801-3

Rims

Description	Color	Order Qty.	Cat. No.	11 Digit Number
-------------	-------	------------	----------	-----------------

Rim, snap-on for

8	Cardiology III (large side)	Black	10/bag	36548	78-8078-2789-0
	Cardiology III (large side)	Gray	10/bag	36549	78-8078-2790-8
9	Cardiology III (small side)	Black	5/bag	36550	78-8078-2791-6
	Cardiology III (small side)	Gray	5/bag	36551	78-8078-2792-4
8	Cardiology II S.E. / Cardiology II	Black	10/bag	36548	78-8078-2789-0
	Cardiology II S.E. / Cardiology II	Gray	10/bag	36549	78-8078-2790-8
8	Classic II S.E. / Classic II	Black	10/bag	36548	78-8078-2789-0
	Classic II S.E. / Classic II	Gray	10/bag	36549	78-8078-2790-8
8	Select	Black	10/bag	36548	78-8078-2789-0
	Select	Gray	10/bag	36549	78-8078-2790-8
8	Lightweight	Black	10/bag	36548	78-8078-2789-0
	Lightweight	Gray	10/bag	36549	78-8078-2790-8
9	Classic II Pediatric	Black	5/bag	36550	78-8078-2791-6
	Classic II Pediatric	Gray	5/bag	36551	78-8078-2792-4
10	Classic II Infant	Black	5/bag	36552	78-8078-2793-2
	Classic II Infant	Gray	5/bag	36553	78-8078-2794-0



Special Procedure Adaptors

Description	Color	Order Qty.	Cat. No.	11 Digit Number
-------------	-------	------------	----------	-----------------

Adaptors for small site or pediatric auscultation

11	Master Cardiology	Black	3/bag	36565	78-8078-2805-4
	Master Cardiology	Gray	3/bag	36566	78-8078-2806-2
12	Cardiology II S.E. / Cardiology II	Black	3/bag	36563	78-8078-2803-9
	Cardiology II S.E. / Cardiology II	Gray	3/bag	36564	78-8078-2804-7

Binaural Assemblies—include tubing, eartubes and eartips


Binaural replacements are available to distributors for use in stethoscope repairs.

Note: Length is measured from tip of eartips to center of chestpiece.

	Description	Color	Length	Part Number
Binaural Assembly for				
1	Master Cardiology	Black	22 in. (55.8cm)	78-8078-6582-5
2	Master Cardiology	Black	27 in. (69cm)	78-8078-6583-3
	Master Cardiology	Navy Blue	27 in. (69cm)	78-8078-6585-8
	Master Cardiology	Burgundy	27 in. (69cm)	78-8078-6586-6
	Master Cardiology	Hunter Green	27 in. (69cm)	78-8078-6587-4
	Master Cardiology	Black Plated, Black	27 in. (69cm)	78-8078-6621-1
1	Cardiology III	Black	22 in. (55.8cm)	78-8078-6582-5
2	Cardiology III	Black	27 in. (69cm)	78-8078-6583-3
	Cardiology III	Navy Blue	27 in. (69cm)	78-8078-6585-8
	Cardiology III	Burgundy	27 in. (69cm)	78-8078-6586-6
	Cardiology III	Hunter Green	27 in. (69cm)	78-8078-6587-4
3	Master Classic II	Black	27 in. (69cm)	78-8062-9927-3
	Master Classic II	Caribbean Blue	27 in. (69cm)	78-8078-8650-8
	Master Classic II	Navy Blue	27 in. (69cm)	78-8062-9744-2
	Master Classic II	Burgundy	27 in. (69cm)	78-8062-9745-9
	Master Classic II	Black Plated, Black	27 in. (69cm)	78-8063-1137-5
	Master Classic II	23 Carat Gold Plated	27 in. (69cm)	78-8063-1142-5
	Master Classic II	Hunter Green	27 in. (69cm)	78-8078-8651-6
	Master Classic II	Purple	27 in. (69cm)	78-8078-5846-5
4	Classic II S.E. / Classic II Pediatric / Classic II Infant	Black	28 in. (71.1cm)	12-2376-8527-1
	Classic II S.E. / Classic II Pediatric / Classic II Infant	Gray	28 in. (71.1cm)	12-2376-8528-9
	Classic II S.E. / Classic II Pediatric / Classic II Infant	Red	28 in. (71.1cm)	12-2376-8532-1
	Classic II S.E. / Classic II Pediatric / Classic II Infant	Caribbean Blue	28 in. (71.1cm)	12-8078-8636-7
	Classic II S.E. / Classic II Pediatric / Classic II Infant	Navy Blue	28 in. (71.1cm)	12-2376-8530-5
	Classic II S.E. / Classic II Pediatric / Classic II Infant	Hunter Green	28 in. (71.1cm)	78-8078-6502-3
	Classic II S.E. / Classic II Pediatric / Classic II Infant	Purple	28 in. (71.1cm)	78-8078-6503-1
	Classic II S.E. / Classic II Pediatric / Classic II Infant	Burgundy	28 in. (71.1cm)	12-2376-8533-9
	Classic II S.E. / Classic II Pediatric / Classic II Infant	Raspberry	28 in. (71.1cm)	78-8078-8635-9

Binaural Assemblies—include tubing, eartubes and eartips (cont.)

Description	Color	Length	Part Number
Binaural Assembly for			
5 Select / Lightweight	Black	28 in. (71.1cm)	78-8078-1206-6
Select / Lightweight	Caribbean Blue	28 in. (71.1cm)	78-8078-8659-9
Select / Lightweight	Raspberry	28 in. (71.1cm)	78-8078-8658-1
Select / Lightweight	Light Blue	28 in. (71.1cm)	12-2376-1252-3
Select / Lightweight	Burgundy	28 in. (71.1cm)	78-8078-4518-1
Select / Lightweight	Royal Blue	28 in. (71.1cm)	78-8078-4520-7
Select / Lightweight	Purple	28 in. (71.1cm)	78-8078-5839-0
Teaching Classic II S.E., 2 headsets (not shown)	Black	40 in. (100cm)	78-8078-0284-4
Teaching Master Classic II, 2 headsets (not shown)	Black	40 in. (100cm)	78-8078-0285-1
Anesthescope (not shown)	Black	48 in. (120cm)	12-2376-3734-8



3M™ Double-Stick Discs

Pre-cut, hypoallergenic and easy to use, 3M Double-Stick Discs are ideal for securing Littmann stethoscope chestpieces to the skin, freeing the user's hands to operate other pieces of equipment.



20

Description	Order Qty.	Cat. No.	Stock No.
Double-Stick Discs	612/case	2181	70-2221-8101-1

Identification Tags

A high quality I.D. tag that is both durable and attractive. The customer has the option to affix either a write-on or engravable name plate (both included in the package).



21

Description	Color	Order Qty.	Cat. No.	Stock No.
3M™ Littmann™ Identification Tag				
Tag fits all sizes of Littmann Stethoscopes	Black	5/case	2170	70-2006-7010-0
	Gray	5/case	2171	70-2006-7011-8
Tag fits all sizes of Littmann Stethoscopes (for customers outside the USA)	Black	5/case	2170X	70-2006-7923-4
	Gray	5/case	2171X	70-2006-7924-2

3M Littmann Stethoscope Warranty and Repair Program

A 3M Littmann stethoscope is a lifetime investment. Beyond the liberal warranty policy, 3M has a lifetime repair and replacement program that is unsurpassed in the industry.

The Littmann Master Cardiology stethoscope is warranted against defects in material and workmanship for a period of seven years from the date of purchase.

The Littmann Cardiology III and Cardiology II S.E. stethoscopes are warranted against defects in material and workmanship for a period of five years from the date of purchase.

The Littmann Master Classic II, Classic II S.E., Classic II Pediatric, Classic II Infant, and all speciality stethoscopes (Littmann Teaching, and Anesthescope stethoscopes), are warranted against any defects in material and workmanship for a period of three years from the date of purchase.

The Littmann Select and Lightweight Stethoscopes are warranted against defects in material and workmanship for a period of two years from the date of purchase. Except in cases of obvious abuse or accidental damage, if a material or manufacturing defect is discovered during the warranty period, repairs will be made without charge upon return of the instrument to 3M.

Return instrument by traceable courier for repair to:

3M Health Care Service Center
 Building 502, Suite 200
 3350 Granada Ave. North
 Oakdale, MN 55128
 (800) 292-6298

To order parts and accessories call or fax:

Phone: (800) 327-5380 — Hospitals
 (800) 592-3979 — Distributors
 Fax: (800) 772-2547

Outside of the USA, please contact your local 3M subsidiary.



Health Care

3M Center, Building 275-4E-01
 St. Paul, MN 55144-1000
 USA
 1 800 228-3957
<http://www.3M.com/Littmann>

3M Canada Inc.

Post Office Box 5757
 London, Ontario N6A 4T1
 Canada
 1 800 563-2921

3M Puerto Rico

Puerto Rico Industrial Park
 P.O. Box 100
 Pueblo Station
 Carolina, Puerto Rico 00986-0100
 1 809 750-3000



40% Pre-consumer waste paper
 10% Post-consumer waste paper

Printed in U.S.A.
 © 3M 2000 70-2008-3362-5 WG