



# PRODUCT SPECIFICATION



**Product:** EMERGENCY BLANKET

**REF:** 1892



## 1 General description

An aluminum coated plastic film packed in a plastic package with printed illustrations. Gold colored on one side, silver on the other. The blanket helps prevent hypothermia, assist in prevention of post-accident shock and protects the injured from wind and moist.

## 2 Intended use

The Emergency Blanket is intended for use in a first aid situation.

## 3 Technical description

### *Product dimensions*

Blanket Unfolded: 210 x 160 cm , folded ~11 x 8 cm.

### *Material*

Film: Vacuum aluminized PET (Polyethylene terephthalate)

Package: PP (polypropylene)

## 4 Composition

### *Consumer package*

<u>Ref#</u>	<u>Product/Material</u>
51618920	Emergency blanket in package

### *Transport package*

<u>Pcs</u>	Consumer packages (pcs):	10
	Dimensions (l x w x h) cm:	21x16x11
	Weight (kg):	0,7

## 5 Labeling

### *On product/ consumer package*

Cederroth address, telephone- and faxnumber, web-/mail-address  
Product name, Reference code  
Lot-number, Expiry date, EAN-code, CE  
Instructions

### *On transport package*

Cederroth address, telephone- and faxnumber, web-/mail-address  
Product name, Reference code  
Lot-number, Expiry date, EAN-code, CE  
Quantity

## 6 Storage conditions

The product should be stored in the original package for maximum protection until time of use.

Shelf life: 5 years from production date

## Authorization: