



CODMAN® Craniotomy Products

The cornerstone of neurosurgery



Codman

Over 30 years of providing quality products!

CODMAN Surgical Patties & Strips

CODMAN Surgical Patties & Strips provide the neurosurgeon with softness, gentleness, protection and absorbency to meet the demands of delicate intracranial and spinal surgery.

Absorbency


CODMAN Patties & Strips absorb more than 5 times their weight in less than a second.

Precision Cut

CODMAN patties and strips precision cut provides smooth edges and helps prevent linting.

COTTONOID® Material

COTTONOID compressed rayon provides the natural softness, absorbency and pliability of cotton. A special process interweaves individual rayon fibers to resist shedding, stretching or linting – whether wet or dry.



CODMAN Competitor

Absorbency comparison dye test



Soft and Gentle

SUTUREWELD®'s Strong bond holds the indicator string to the patty so the surface which lies next to tissue is completely smooth and free of stitches.

Pliable and Porous

COTTONOID material offers the advantage of being porous so CODMAN Patties & Strips easily conform to the intricate contours of delicate tissue. Yet, you can safely suction through them without aspirating the tissue or the patty or strip.

X-Ray Detectable

Marking is where it should be – on the patty or strip.

Unique Packaging

A transparent, sterile envelope shows the size and shape of patties and strips, and the Count Verification Card instantly shows how many patties you started with and how many were used.

To Order Call 1-800-225-0460

www.codman.com

See back page for sizes and catalog numbers.

BICOL® Collagen Sponge – Proven Quality

BICOL Collagen sponge is a highly useful adjunct to neurosurgery as a moistening agent and as a protective device. The inherent composition of the material provides exceptional wet strength while maintaining its soft, pliable and resilient nature.



80-1460 10 per box

Moistening

As a moistening agent, the BICOL Collagen Sponge is highly efficient during prolonged periods of exposure. The material can rapidly absorb and maintain up to 60 times its weight in aqueous fluid, in comparison to the standard gelatin sponge, UPS, which absorbs about 30 times its weight. Additionally, its microfibrillar structure sustains its pliability and prevents moistened material from crumbling or breaking apart while handling.

Protection

As a protective device the resilient nature of the moistened material effectively buffers the pressure of a retractor blade. After its rapid initial wetting, there is little change in its overall volume and no change with the passage of time following application.

Application

Application and manipulation of the wetted material is facilitated by its surface textures. One surface is very smooth so that once moistened it slides over the brain and readily conforms to the area to be protected. The opposite surface has sufficient grain to provide friction to a retractor blade for manipulation.

CRANIOPLASTIC™ Kit (43-1050)

CRANIOPLASTIC is a resinous material for repairing cranial defects.

Two complete sets consist of two bags sterile powder (net weight 30 grams each), two vials sterile liquid (net contents 17 milliliters each) and two sterile matrix bags.

ANEUROPLASTIC™ Kit (43-1060)

Aneuroplastic is a fast-setting resinous material which is used to encapsulate cerebral aneurysms.

Two complete sets consist of two bags sterile powder (net weight 30 grams each), two ampules sterile liquid (net contents 17 milliliters each), and two sterile matrix bags.



To Order Call 1-800-225-0460
www.codman.com

See back page for sizes and catalog numbers.

CODMAN Raney Scalp Clips & Applier

A complete system for the neurosurgeon

CODMAN Disposable Scalp Clips:

- Provides for hemostasis of scalp flap edges

CODMAN Scalp Clip Forceps

Single-use clips used with:


- RANEY Scalp Clip Applying Forceps (20-5115)
- LEROY (LEROY-RANEY) Scalp Clip Applying Forceps (20-5117)
- Infant LEROY Scalp Clip Applying Forceps (20-5120)

The CODMAN Scalp Clip Applier is a sterile, single-use system designed to apply scalp clips to the scalp flap for hemostasis during a craniotomy procedure.

- Comes with pre-loaded scalp clips
- Permits rapid application of scalp clips
- Provides adequate hemostatis



LEROY-RANEY Scalp Clips 20-1034 (Used with the CODMAN Disposable Scalp Clip Applier)



Squeeze the scalp clip applier to load scalp clip into position

Applier is now engaged to open clip to position on scalp flap and is quickly and easily disengaged.

Applier allows maneuverability for all clip positions.

Scalp clips are removed with scalp clip removal forceps.

CODMAN Surgical Marker (43-1030, 36/case) CODMAN Microsurgical Marker (43-1032, 36/case)

CODMAN now offers two surgical marking pens to meet the need for both a bold and fine line marker.

The CODMAN Surgical Marker delivers a traditional bold line while the new CODMAN Microsurgical Marker offers a particularly fine line for plastic and other delicate procedures.

The tip on the CODMAN Microsurgical Marker is slightly flexible.



Important Features

- Both pens are completely filled with pre-resistant gentian violet ink to reduce the chance that they will dry out in storage.
- Individually packaged, sterile, in see-through peel pouches.
- Dispenser boxes make storage convenient.
- Each pen is packaged with a 15cm/6" disposable ruler.

To Order Call 1-800-225-0460
www.codman.com

See back page for sizes and catalog numbers.

Sterile! Sharp!

CODMAN Disposable Perforator Family

Sterile, sharp, disposable perforators.

Designed for use in perforating the cranium.

100% Quality Testing: All units are pretested to meet specifications.

Available in three convenient sizes; 9mm, 11mm, and 14mm

Also available in Craniotomy Kits

Perforators are color-coded for easy identification



CODMAN Craniotomy Kit

CODMAN brings you a complete craniotomy system in a convenient three-piece kit. Included for each size are a CODMAN Disposable Perforator, CODMAN Cranio-blade and CODMAN Wire Pass Drill with the following features:

CODMAN Disposable Perforators*

- Single use, sterile
- Color-coded sizes
- No maintenance required
- Sharp

CODMAN Cranio-blade

- Four-fluted channel
- Single use
- Sterile

CODMAN Wire Pass Drill

- Single use
- Sterile



26-1246



26-1247



*Available in 9mm, 11mm, or 14mm

To Order Call 1-800-225-0460

www.codman.com

See back page for sizes and catalog numbers.

Ordering Information

CODMAN Surgical Patties

Sterile, x-ray detectable, 10 patties per package, 20 packages per dispenser box

Order #	Description
80-1396	5/16" Round (.79cm)
80-1397	1/4" x 1-1/2" (64cm x 3.81cm)
80-1398	1/4" x 3" (64cm x 7.62cm)
80-1399	1/4" x 1/4" (64cm x 64cm)
80-1400	1/2" x 1/2" (1.27cm x 1.27cm)
80-1401	3/4" x 3/4" (1.91cm x 1.91cm)
80-1402	1/2" x 1" (1.27cm x 2.54cm)
80-1403	1" x 1" (2.54cm x 2.54cm)
80-1404	1/2" x 1-1/2" (1.27cm x 3.91cm)
80-1406	1/2" x 2" (1.27cm x 5.06cm)
80-1407	1/2" x 3" (1.27cm x 7.62cm)
80-1408	1" x 3" (2.54cm x 7.62cm)
80-1409	3" x 3" (7.62cm x 7.62cm)

CODMAN Surgical Strips

Sterile, x-ray detectable, 10 strips per package, 20 packages per dispenser box

80-1449	1/8" x 6" (3mm x 152mm)
80-1450	1/4" x 6" (6mm x 152mm)
80-1451	1/2" x 6" (13mm x 152mm)
80-1452	3/4" x 6" (19mm x 152mm)
80-1453	1" x 6" (26mm x 152mm)
80-1454	1-1/2" x 6" (38mm x 152mm)
80-1455	2" x 6" (51mm x 152mm)
80-1456	3" x 6" (76mm x 152mm)
80-1457	3-1/2" x 6" (89mm x 152mm)

CODMAN Disposable Perforator with Hudson son

Individually sterilized, one per package

26-1221	CODMAN Disposable Perforator, 14mm
26-1222	CODMAN Disposable Perforator, 11mm
26-1223	CODMAN Disposable Perforator, 9mm

CODMAN Craniotomy Kit Sterile, disposable

Consists of one each: CODMAN Disposable Perforator, CODMAN Cranio-blade, and CODMAN Wire Pass Drill

26-1230	CODMAN Craniotomy Kit, 14mm
26-1231	CODMAN Craniotomy Kit, 11mm
26-1229	CODMAN Craniotomy Kit, 9mm

26-1246	CODMAN Cranio-blade Sterile, disposable, five per package
26-1247	CODMAN Wire Pass Drill Sterile, disposable, five per package

80-1460	BICOL Collagen Sponge BICOL sponges are 3" (76mm) x 4" (102mm) and packaged individually, sterile, ten packages per box
---------	---

20-5115	RANEY Scalp Clip Applying Forceps Use with 20-1037 or 20-1025 RANEY Scalp Clips Length 6-1/4" (157mm)
20-5117	LEROY Scalp Clip Applying Forceps Adult size for use with 20-1032 or 20-1034 LEROY Scalp Clips Length 7-1/4" (184mm)
20-5120	Infant LEROY Scalp Clip Applying Forceps For use with 20-1035 Infant LEROY Scalp Clips Length 7-1/4" (184mm)
20-1025	RANEY Scalp Clips Stainless steel, twelve clips per package Use with 20-5115 RANEY Scalp Clip Applying Forceps
20-1032	LEROY Scalp Clips Sterile, disposable, plastic Ten clips per package, twenty packages per case Use with 20-5117 LEROY Scalp Clip Applying Forceps
20-1034	LEROY-RANEY Scalp Clips Sterile, disposable, plastic Ten clips per package, twenty packages per case Use with 20-5117 LEROY Scalp Clip Applying Forceps
20-1035	Infant LEROY Scalp Clips Sterile, disposable, plastic Ten clips per package, twenty packages per case Use with 20-5120 Infant LEROY Scalp Clip Applying Forceps
20-1037	RANEY Scalp Clips Sterile, disposable, plastic Ten clips per package, twenty packages per case Use with 20-5115 RANEY Scalp Clip Applying Forceps
20-1038	CODMAN Disposable Scalp Clip Applier Sterile, plastic Ten packages per case
20-1039	CODMAN Scalp Clip Removable Forceps Use with 20-1038 CODMAN Disposable Scalp Clip Applier

Order #	Description
43-1006	SPENCE Cranioplastic Roller

43-1022	Titanium Mesh Sterile, titanium mesh is indicated for use with methyl methacrylate formulated for cranioplasty Size 6" x 6" (15cm x 15cm) Wire diameter 0.010" (31g) Mesh per inch 24 x 24
---------	---

43-1050	CRANIOPLASTIC Kit Two complete sets consist of two bags sterile powder (net weight 30 grams each), two vials sterile liquid (net contents 17 milliliters each), and two sterile matrix bags
---------	---

43-1060	ANEUROPLASTIC Kit Two complete sets consist of two bags sterile powder (net weight 30 grams each), two ampules sterile liquid (net content 17 milliliters each), and two sterile matrix bags
---------	--

80-1163	MALIS Integrated Irrigation Tubing & Bipolar Cord Set Sterile, disposable, ten per case
---------	---


80-9102	MALIS Integrated irrigation Tubing & Bipolar Cord Set Length 10 ft (for Module 110 Irrigator)
---------	--

80-9103	Disposable Irrigation Tube Set
---------	--------------------------------

30-1539	Bipolar Cords Sterile, disposable, ten per case, 12 ft
---------	---

General Disposables

19-4137	Towel Clamp, plastic, disposable, 6/package, 24 per case
43-1030	CODMAN Surgical Marker, 36 per case
43-1032	CODMAN Microsurgical Marker, 36 per case
63-4031	Vein Stripper, 12 per case



Codman
a Johnson & Johnson company

325 Paramount Drive
Raynham, MA 02767 USA

For more information, contact your CODMAN Sales Representative
For product information, call 800-225-0460

www.codman.com