

## Extension set

20350E-0006



- Needle-Free Connectors
- Set Length 43 cm
- Priming Volume 5 ml
- Bore 2.7 mm
- Vented Filter
- 1 SmartSite® Valve Y Port
- 1 Pinch Clamp(s)
- 0.2 µm Inline Filter
- 1 Male Luer(s)
- Rotating Male Luer(s)
- Low Sorbing
- Polyethylene lined PVC
- Natural rubber latex is not part of the material formulation
- DEHP is not part of the material formulation
- Quantity 100

